




ENDOSCREENER COMPUTER AIDED POLYP DETECTION SYSTEM

ARTIFICIAL INTELLIGENCE CAD_e



CONSISTENT POLYP DETECTION – ALL THE TIME

Human Central Vision Has A Tight Focus Of Around 13°¹
Computer Vision Does Not Have That Limitation
EndoScreener[®] Focuses On 100% Of The Screen - 100% Of The Time

Always On – Always Focused – Doesn't Get Tired

EndoScreener[®] software utilizes deep learning artificial intelligence (A.I.) and Computer Aided Polyp Detection (CAD_e) systems to efficiently analyze and identify local abnormalities in the lining of the colon, significantly improving the detection and precision capabilities of a standard colonoscopy.

EndoScreener, a **software as a medical device (SaMD)** product, is compatible with most mainstream endoscope systems.

TECHNOLOGY IS A GREAT LEVELER

Introducing EndoScreener - A Technologically Advanced AI Polyp Detection System

Conventional AI systems employ an **Object Detection Approach**, that is focused on a lesion's general appearance.

EndoScreener utilizes advanced mathematical modeling and **Pixel Level, Image Segmentation** to produce an algorithm focused on **Local Features**

Pixel level awareness of local features can improve an AI polyp detection system's ability to be **more sensitive to difficult lesions** while producing **fewer false positives**²



EndoScreener Algorithm Segments and Identifies **All** Pixels In A Frame That May Be Associated with a Polyp Before Drawing The "Alert Box"

In clinical use, a local-feature oriented algorithm, like EndoScreener, considers details that are less obvious to the human eye and may be more effective in difficult cases, such as:

- Sessile serrated lesions that have fewer macro features
- Obscured / partially appearing lesions
- Flat lesions

STRIKING A BALANCE BETWEEN USEFUL AND DISTRACTING

How well does AI Find the target?

Sensitivity is a measure of a systems ability to identify the presence of a polyp in the visual field

Frame based sensitivity is the most important metric for an AI polyp detection system



	GI Genius™	EndoScreener
Per Frame Sensitivity	47.73% ³	91.64% ²

ENDOSCREENER'S SUPERIOR SENSITIVITY MAY HELP TO:

- Find difficult to detect flat, subtle, and laterally spreading lesions
- Alert the user to partially appearing polyps
- Alert the user sooner when the polyp enters the visual field

FEWER FALSE POSITIVES

- EndoScreener Focuses On Polyps
Not Non-polyp Objects Like Folds, Feces & Bubbles

	GI Genius™	EndoScreener
False Alarms Per Colonoscopy	101 ³	26 ⁴

PROVEN CLINICAL PERFORMANCE – SUPPORTS HIGH QUALITY OUTCOMES

- More Than 5,000 Patients In Six Rigorous, Randomized Controlled Trials (RCTs)
- Two Tandem Colonoscopy Randomized Controlled Trials
- Double Blind Randomized Controlled Trial


In A Back-To-Back, Tandem Study, EndoScreener Increased Adenoma Detection Rate (ADR) by 8.8%.
A 43.4% Relative Increase Compared With Traditional White-Light Colonoscopy⁵

A 1% increase in ADR results in a 3% decrease in the risk of interval cancer and a 5% decrease in the risk of a fatal interval colorectal cancer⁶


“Sessile Serrated Lesion (SSL) miss rate was significantly reduced in our US Trial.

The world’s 1st Randomized Controlled Trial That Demonstrated Statistically Measurable Improvement In SSL Metrics⁷


ADR Improvement⁵



Improved Lesion Detection⁵




Sessile Serrated Lesion Miss Rate⁷



EndoScreener use resulted in 35% fewer missed Sessile Serrated Lesions in physicians with a high baseline ADR⁷

Adenoma Miss Rate⁷



EndoScreener use resulted in 11% fewer missed adenomas in physicians with a high baseline ADR (44%)⁷


COLONOSCOPY COMBINED WITH AI MAY SIGNIFICANTLY INCREASE FOLLOW UP VISITS

The use of artificial intelligence during colonoscopy increased the proportion of patients requiring intensive colonoscopy surveillance by approximately 35% in the United States (absolute increases of 2.9%)⁸


CHALLENGING CASES DISCOVERED BY ENDOSCREENER

EndoScreener's Segmentation Technology Focuses On **Local Features** That Are Less Obvious To The Human Eye.


This Local Feature Approach May Be More Effective In Finding Partially Obscured Lesions While Producing Fewer False Positives.




Conventional



Conventional




Conventional




Flat & Isochromatic



Partly behind colon folds




Insufficient air inflation



Unclear boundary



Overexposure



Partial occlusion – liquid



Blurred lens



Insufficient light



Partial occlusion – feces

References:

- ¹Yonis O, et al. A Hazard Detection and Tracking System for People with Peripheral Vision Loss using Smart Glasses and Augmented Reality. (IJACSA) International Journal of Advanced Computer Science and Applications, Vol. 10, No. 2, 2019 DOI: 10.14569/IJACSA.2019.0100201
 - ²Wang P, Xiao X, Glissen Brown JR, et al. Development and validation of a deep-learning algorithm for the detection of polyps during colonoscopy. Nat Biomed Eng. 2018;2(10):741-748. doi:10.1038/s41551-018-0301-3
 - ³Brand, M et al. Frame-by-Frame Analysis of a Commercially Available Artificial Intelligence Polyp Detection System in Full-Length Colonoscopies, Digestion. Published online: June 29, 2022. DOI: 10.1159/000525345
 - ⁴Holzswanger EA, Bilal M, Glissen Brown JR, Singh S, Becq A, Ernest-Suarez K, Berzin TM. Benchmarking definitions of false-positive alerts during computer-aided polyp detection in colonoscopy. Endoscopy. 2021 Sep;53(9):937-940. doi: 10.1055/a-1302-2942. Epub 2021 Jan 18. PMID: 33137833; PMCID: PMC8386281.
 - ⁵Wang P, Berzin TM, Glissen Brown JR, Bharadwaj S, Becq A, Xiao X, Liu P, Li L, Song Y, Zhang D, Li Y, Xu G, Tu M, Liu X. Real-time automatic detection system increases colonoscopic polyp and adenoma detection rates: a prospective randomised controlled study. Gut. 2019 Oct;68(10):1813-1819. doi: 10.1136/gutjnl-2018-317500. Epub 2019 Feb 27. PMID: 30814121; PMCID: PMC6839720.
 - ⁶Corley DA, Jensen CD, Marks AR, Zhao WK, Lee JK, Doubeni CA, Zauber AG, de Boer J, Fireman BH, Schottinger JE, Quinn VP, Ghai NR, Levin TR, Queensberry CP. Adenoma detection rate and risk of colorectal cancer and death. N Engl J Med. 2014 Apr 3;370(14):1298-306. doi: 10.1056/NEJMoa1309086. PMID: 24693890; PMCID: PMC4036494.
 - ⁷Glissen Brown JR, Mansour NM, Wang P, et al. Deep Learning Computer-aided Polyp Detection Reduces Adenoma Miss Rate: A United States Multi-center Randomized Tandem Colonoscopy Study (CADeT-CS Trial). Clin Gastroenterol Hepatol. 2022;20(7):1499-1507.e4. doi:10.1016/j.cgh.2021.09.009
 - ⁸Mori, Y et al. Impact of Artificial Intelligence on Colonoscopy Surveillance After Polyp Removal: A Pooled Analysis of Randomized Trials. Clinical Gastroenterology & Hepatology Published: August 26, 2022, DOI: <https://doi.org/10.1016/j.cgh.2022.08.022>
- "GI Genius™" is a trademark of Medtronic

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTEndoscopy.com
customer@mtendoscopy.com

MPW-WIPD01-B Rev 2