

SURECLIP

ACCURATE & EFFECTIVE HEMOSTASIS




SureClip[®] is designed to provide precise clip placement and improved visualization.

Unrestricted opening and closing of the jaw prior to deployment

Reliable 360° bidirectional rotation, even in tortuous anatomy


Shorter stem makes the clip less obtrusive and improves visualization of the target area

Available in four different sizes (8 mm, 11 mm, 16 mm, and 17 mm)



GREATER CLOSURE/ COMPRESSION FORCE¹


Testing was completed to determine the compression force between the clip arms and overall locking force. The closure/compression force evaluation confirmed that both SureClip and SureClipPlus have greater average clamping strength than that of Boston Scientific Resolution[®].



EXCELLENT RETENTION STRENGTH²

When comparing the SureClip^{PLUS} to other competitive hemostasis clips, testing showed that SureClip^{PLUS} has greater retention strength on average than Boston Scientific Resolution[®]. Additionally, SureClip demonstrated similar strength on average compared to Resolution[®].

The SureClip's rotatability is smooth, with a handle that allows controlled movements with intuitive functionality.



“The SureClip is, in my opinion, clearly better than the competitions' clip because of its **accuracy, reliability, and functionality.**”

- VK Sharma, MD, FASGE, FACG, AGAF

NICKEL (NI) SENSITIVITY³


Several metals may be used in the manufacture of hemostasis clips and most have varying degrees of nickel content. SureClip is made from medical-grade stainless steel, which contains some nickel.

MATERIAL / APPROXIMATE % OF NICKEL

Stainless Steel	Elgiloy	Nitinol
6-8%	14-16%	>50%

When nickel is present, the surface area of the clip and its potential contact with patient may also be relevant. In testing, SureClip was found to have the smallest surface area when compared with 3 competitive hemostasis clips.

ESTIMATED TOTAL SURFACE AREA



SPECIFICATIONS

SURECLIP^{MINI} LOWEST PROFILE STERILE

Order Number	Henry Schein Item Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size(mm)	Package Units
RC30415	132-5180	8	Max 2.6	235	2.8	2/Box
RC30411	132-5187	8	Max 2.6	235	2.8	10/Box

SURECLIP SHORTEST STEM STERILE

RC30445	132-5724	11	Max 2.6	235	2.8	2/Box
RC30441	132-5723	11	Max 2.6	235	2.8	10/Box

SURECLIP^{PLUS} GREATEST CLOSURE/COMPRESSION FORCE STERILE

RC30385	128-5657	16	Max 2.6	235	2.8	2/Box
RC30381	128-5655	16	Max 2.6	235	2.8	10/Box

SURECLIP^{MAX} WIDEST JAW OPENING* STERILE

RC30545	136-6309	17	Max 2.6	235	2.8	2/Box
RC30541	136-6310	17	Max 2.6	235	2.8	10/Box

*Jaws will open to the maximum reach when placed in contact with tissue.

INDICATIONS FOR USE

The SureClip Repositionable Hemostasis Clip is indicated for Endoscopic clip placement within the Gastrointestinal tract of adult patients only via a straight or side-viewing flexible endoscope for the purpose of :

- (1) Endoscopic marking;
- (2) Hemostasis for
 - (a) Mucosal / sub-mucosal defects < 3 cm,
 - (b) Bleeding ulcers,
 - (c) Polyps < 1.5 cm in diameter,
 - (d) Diverticula in the colon,
 - (e) Arteries < 2 mm
 - (f) Prophylactic clipping to reduce the risk of delayed bleeding post lesion resection;
- (3) As a supplementary method, closure of GI tract luminal perforations <20 mm that can be treated conservatively.

MR CONDITIONAL

Non-clinical testing has demonstrated the SureClip Repositionable Hemostasis Clip is MR Conditional according to ASTM F2503.

¹ Bench testing was performed (n=12) per device. Test report: MTN-ROCC-TR/JP-004.

² Bench testing was performed (n=20) per device. Test report: MTN-ROCC-TR-CM-01 Rev:A2.

³ Independent testing was conducted at the Modern Analysis Center, Nanjing University.

All trademarks are the property of their respective owners.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTEndoscopy.com
customercare@mtendoscopy.com

MPW30000R Rev.5