Cold Snare

Instructions for Use

Asa fría

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Alça a Frio

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Instruction page for Electronic instructions for use website:

1. The eIFU is only applicable for Brazil.
2. Reference the product labeling and choose the correct version of eIFU.

To access the instructions for use, please go to <https://eifu.micro-tech.com.cn> by entering the product-specific IFU code CS05. For a paper copy, provided in 7 days at no cost, please call (+1) 877-552-4027.

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IMPORTANT INFORMATION

NE

Caution: Federal law restricts this device to sale by or on the order of a physician. Read all instructions carefully before use. They contain essential information on using this device safely and effectively. Keep these instructions in a safe, accessible location, as you may need to refer to them again. If you have any questions or comments about any information in these instructions, please contact Micro-Tech.

INTENDED USE

Cold Snare is intended to be used without electro-cautery in the gastrointestinal tract for the endoscopic resection of diminutive polyps($\leq 9\text{mm}$):

1. Single or multiple polyps in digestive tract, the diameter should generally being $\leq 9\text{mm}$, the number of the multiple polyps being ≤ 30 ;
2. Pathologic histology confined that polyp is non- infiltration;
3. Early cancer confined to a mucosa and the polyp(s) suitable for endoscopic removal

CONTRAINDICATIONS

Contraindications include, but may not be limited to:

1. People who have poor physical fitness, and a serious heart or/and lung disease, and cannot tolerate to endoscopy and endoscopic treatment.
2. People who have a bleeding tendency, an extended bleeding and coagulation time, or thrombocytopenia or prolonged prothrombin, which cannot be corrected by treatment.
3. People who have a polyp / polyps with a too large base, the base of a polyp $> 9\text{mm}$.
4. People who have a polypoid cancer which has infiltrated a tissue or organ and has deteriorated.
5. People who have diabetes, regardless of whether whose blood sugar is normal should be classified as relative contraindications (poor healing capacity).
6. Uncooperative patients or family members.

POTENTIAL COMPLICATIONS

1. Bleeding, perforation, infection, septicemia, etc.
2. Respiratory depression or apnea, arrhythmia or cardiac arrest.
3. Potential complications that are not known or observed.

WARNINGS


1. **The product is intended for single use only! DO NOT re-use, re-sterilize, and/or reprocess.** Re-use re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed.
2. Do not use this instrument for any purpose other than its intended use.
3. The product is only intended for adult populations.
4. This device is not made with natural rubber latex.
5. Patients should be informed of the potential risks and complications, which may lead to injury, illness or death of the patient.
6. **The instrument is intended for use under the direct supervision of a suitably trained physician only.** A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
7. Confirm that the endoscopy view is clear before use. Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. Damage to the endoscope and/or the instrument may also occur.
8. Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view. If you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or instrument.
9. The operator of the snare must be a physicians with sufficient training.
10. If polyps cannot be removed, and don't force operation, Please changing cold snare or using other methods.

11. For multiple polyps patients, after repeated procedures, please observe whether loop is distorted or not, if distortion, please use another one.
12. For less than 3mm polyps in the digestive tract is not recommended.

- 【 Product Name 】 Cold Snare
- 【 Packaging 】 Packed in Pouch
- 【 Production Date 】 See packaging
- 【 Sterilization 】 Sterilized by EO (ethylene oxide) gas
- 【 Shelf Life 】 5 years
- 【 Compatible Working Channel 】 $\geq \phi 2.8 \text{ mm}$

STRUCTURE

Cold Snare mainly includes a snare loop, wire cable, outer tube and handle. (Fig. 1).

1. Snare Loop 2. Wire Cable 3. Outer Tube 4. Handle Assembly

Fig.1 Schematic diagram of cold snare

Note: This is the sketch map of one kind of loop, the other loop shapes are not shown here.

PREPARATION

1. Reference the product label and choose the appropriate device.
2. Contents supplied STERILE.
3. Inspect the package before use for any damage. Do not use if package is damaged. Open the package carefully after verifying the shelf life.
4. Carefully remove the device from its packaging and uncoil it. Do not use excessive force as this may damage the device and affect performance.
5. Before use, please check the insertion patients body parts, ensure that no sharp edge.
6. If this device shows any signs of damage, do not use. Do not attempt to repair a nonfunctional or damaged device.

EN

INSTRUCTIONS FOR USE

1. Cold Snare is compatible with an endoscope channel of 2.8mm or larger.
2. Retract the handle in order to make sure the snare loop is inside the outer tube (as shown in Fig.2). And insert the snare into the working channel of the endoscope. Via the endoscope view, make sure that the tube projects out of the distal end of the channel (as shown in Fig. 3).


Fig. 2


Fig. 3

3. Extend the snare loop, trap the polyps gently (as shown in Fig. 4).


Fig. 4

4. Pull back the slider to removing polyps.
5. Retract the handle in order to make sure the snare loop is inside the outer tube , and pull out cold snare Repeat the procedure as necessary.
6. Upon completion of the procedure, remove the instrument from the endoscope.

STORAGE

The product should be stored in a cool, dry, clean, well-ventilated, non-corrosive gas environment.

Do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

The product shelf life is 5 years.

PRODUCT DISPOSAL

After use, this product may be a potential biohazard. Handle and dispose of in accordance with hospital, local and administrative laws and regulations.

Limited Warranty and Disclaimers

1. Limited Warranty to Buyer. Micro-Tech USA warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored

and used in accordance with the instructions for storage and use provided by Micro-Tech USA and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech USA's literature are meant to generally describe the products and do not constitute any express warranties. In the event that Micro-Tech USA gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech USA's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech USA in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

2. Disclaimer and Release. **THE WARRANTIES, OBLIGATIONS, AND LIABILITIES OF Micro-Tech USA AS SET FORTH HEREIN ARE EXCLUSIVE. BUYER HEREBY WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY BUYER, INCLUDING, BUT NOT LIMITED TO: (1) ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND (2) ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING, OR USAGE OF TRADE.**

3. Implied Warranties. The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech USA to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.

4. Limitation of Liability. **EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, Micro-Tech USA'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO: (a) THE REPLACEMENT OF THE PRODUCTS OR THE RE-SUPPLY OF EQUIVALENT PRODUCTS; (b) THE REPAIR OF THE PRODUCTS OR PAYMENT OF THE COST OF REPAIRING THE PRODUCTS; or (c) PAYMENT OF THE COST OF REPLACING THE PRODUCTS OR ACQUIRING EQUIVALENT PRODUCTS. MICRO-TECH USA SHALL HAVE NO OBLIGATION OR LIABILITY, WHETHER ARISING IN CONTRACT (INCLUDING WARRANTY), TORT (INCLUDING ACTIVE, PASSIVE, OR IMPUTED NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY) OR OTHERWISE, FOR ANY SPECIAL, CONSEQUENTIAL, PUNITIVE, INCIDENTAL, OR INDIRECT DAMAGES, OR FOR LOSS OF USE, LOSS OF REVENUE, LOSS OF BUSINESS, LOST PROFIT, OR OTHER FINANCIAL LOSS ARISING OUT OF OR IN CONNECTION WITH ANY PRODUCT OR OTHER GOODS OR SERVICES FURNISHED BY MICRO-TECH USA, EVEN IF MICRO-TECH USA WAS AWARE OF THE POSSIBILITY OF SUCH DAMAGES OR LOSS.**

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech USA Inc.