MICRO-TECH™
ENDOSCOPY

**Sphincterotome / short-wire
compatible**
Instruction for Use

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STERILE EO

WARNING

Federal law restricts this device to sale by or on the order of a physician. Contents supplied sterile using ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call manufacturer.

For single use only. Do not re-use, reprocess or re-sterilize. Re-use, reprocessing or re-sterilization may compromise the structural integrity of the device and /or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of produce and packaging in accordance with hospital, administrative and/or local government policy.

The device is intended for adults.

INTENDED USE

The Device is intended to be used in the selective cannulation of the biliary ducts and the transendoscopic sphincterotomy of the Papilla of Vater and/or the Sphincter of Oddi. The device can also be used to inject contrast medium.

CONTRAINDICATIONS

Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic sphincterotomy (ES), but may not be limited to:

1. People who have obstruction of upper and lower digestive tract, the products may not reach the duodenum drop location.
2. Dysfunction of coagulation.
3. Patient has suffered a severe allergic reaction to contrast medium.
4. People who have poor physical fitness, and a serious heart or/and lung disease, and cannot tolerate to endoscopy and endoscopic treatment.
5. Acute pancreatitis and acute biliary tract infection.

INSTRUCTION MANUAL

This instruction manual contains essential information on using the instrument safely and effectively.

Please review this manual thoroughly before using the device. If you have any questions or comments about any information in this manual, please contact manufacturer.

ADVERSE EVENTS

Potential complications include, but are not limited to:

- ◆ Pancreatitis
- ◆ Perforation
- ◆ Hemorrhage
- ◆ Hematoma
- ◆ Cholangitis
- ◆ Septicemia/infection
- ◆ Allergic reaction to contrast medium
- ◆ Fulguration
- ◆ Burns
- ◆ Stimulation
- ◆ Cardiac arrhythmias

WARNINGS

This device is not intended to be used in the presence of Flammable liquids, in an Oxygen Enriched Atmosphere or in the presence of explosive gases. Any electrosurgical device constitutes a potential electrical hazard to the patient and/or operator.

Avoid High Frequency output settings where maximum output voltage exceed rated accessory voltage.

【Product Name】 Sphincterotome / short-wire compatible

【Packaging】 Packed in Pouch

【Production Date】 See packaging

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Shelf Life】 3 years

This device is not made with natural rubber latex.

Applicable Electrosurgical Unit

High-frequency electrosurgical unit which is legal listing in America is recommended such as ERBE Electrosurgery;
Recommended power of High-frequency electrosurgical unit: 30 watts to 70 watts

Rated Accessory Voltage of High-frequency electrosurgical unit:

CUT mode: 1200Vp (2400 Vp-p)

DO NOT use higher repeated peak voltage than the Cut-1200Vp (2400Vp-p)

Note: The other generators with higher rated voltage than the Cut-1200Vp (2400Vp-p) is not recommended.

SPECIFICATION

NO.	REF	Tip Length Tolerance: + 2/0 (mm)	Cutting Wire Length Tolerance: ±2 (mm)	With Insulation Coat	Preloaded Guidewire Length (mm)	Cutting Wire Type	Tip type	Compatible Endoscopes Minimum Insertion Channel Width (mm)
1	RDSP-0520-111	5	20	NO	N/A	Single	Cone	≥3.2
2	RDSP-0520-211	5	20	YES	N/A	Single	Cone	≥3.2
3	RDSP-0520-121	5	20	NO	2600	Single	Cone	≥3.2
4	RDSP-0520-221	5	20	YES	2600	Single	Cone	≥3.2
5	RDSP-0520-131	5	20	NO	4500	Single	Cone	≥3.2
6	RDSP-0520-231	5	20	YES	4500	Single	Cone	≥3.2
7	RDSP-0525-111	5	25	NO	N/A	Single	Cone	≥3.2
8	RDSP-0525-211	5	25	YES	N/A	Single	Cone	≥3.2
9	RDSP-0525-121	5	25	NO	2600	Single	Cone	≥3.2
10	RDSP-0525-221	5	25	YES	2600	Single	Cone	≥3.2
11	RDSP-0525-131	5	25	NO	4500	Single	Cone	≥3.2
12	RDSP-0525-231	5	25	YES	4500	Single	Cone	≥3.2
13	RDSP-0530-111	5	30	NO	N/A	Single	Cone	≥3.2
14	RDSP-0530-211	5	30	YES	N/A	Single	Cone	≥3.2
15	RDSP-0530-121	5	30	NO	2600	Single	Cone	≥3.2
16	RDSP-0530-221	5	30	YES	2600	Single	Cone	≥3.2
17	RDSP-0530-131	5	30	NO	4500	Single	Cone	≥3.2

18	RDSP-0530-231	5	30	YES	4500	Single	Cone	≥3.2
19	RDSP-0520-112	5	20	NO	N/A	Weave	Cone	≥3.2
20	RDSP-0520-212	5	20	YES	N/A	Weave	Cone	≥3.2
21	RDSP-0520-122	5	20	NO	2600	Weave	Cone	≥3.2
22	RDSP-0520-222	5	20	YES	2600	Weave	Cone	≥3.2
23	RDSP-0520-132	5	20	NO	4500	Weave	Cone	≥3.2
24	RDSP-0520-232	5	20	YES	4500	Weave	Cone	≥3.2
25	RDSP-0525-112	5	25	NO	N/A	Weave	Cone	≥3.2
26	RDSP-0525-212	5	25	YES	N/A	Weave	Cone	≥3.2
27	RDSP-0525-122	5	25	NO	2600	Weave	Cone	≥3.2
28	RDSP-0525-222	5	25	YES	2600	Weave	Cone	≥3.2
29	RDSP-0525-132	5	25	NO	4500	Weave	Cone	≥3.2
30	RDSP-0525-232	5	25	YES	4500	Weave	Cone	≥3.2
31	RDSP-0530-112	5	30	NO	N/A	Weave	Cone	≥3.2
32	RDSP-0530-212	5	30	YES	N/A	Weave	Cone	≥3.2
33	RDSP-0530-122	5	30	NO	2600	Weave	Cone	≥3.2
34	RDSP-0530-222	5	30	YES	2600	Weave	Cone	≥3.2
35	RDSP-0530-132	5	30	NO	4500	Weave	Cone	≥3.2
36	RDSP-0530-232	5	30	YES	4500	Weave	Cone	≥3.2

PRECAUTIONS

The Sphincterotome shall only be used by physicians or under the supervision of physicians trained in ERCP or ES. A thorough understanding of the technical principles, Clinical applications and risks associated with ERCP/ES is necessary before using this device.

Safe and effective electrosurgery is dependent not only on equipment design, but also, to a large extent, on factors under the control of the operator. It is important that the following be read, understood, and followed in order to enhance safety and effectiveness:

1. Monitor sphincterotomy under fluoroscopy.
2. Avoid bowing of the sphincterotome while it is in a coiled position.
3. Any use of this device, other than those listed in this instruction, is not recommended.
4. Recommend HF (high frequency) Generator model ERBE electrosurgical.
5. Inspect the device before use. If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit

proper working condition, do not use. Please notify MT for return authorization.

6. The Rated Accessory Voltage for this device is 1200Vp (2400Vp-p).
7. This device is not intended to be used in the presence of Flammable liquids, in an Oxygen Enriched Atmosphere or in the presence of explosive gases. Any electrosurgical device constitutes a potential electrical hazard to the patient and/or operator.
8. This device complies with the EMC standard (IEC 60601-1-2:2014) for medical electrical equipment; When connecting to an instrument that complies with the EMC standard for medical electrical equipment, the whole system complies with the EMC standard for medical electrical equipment.
9. The output power selected should be as low as possible for the intended purpose. Certain devices or accessories may present an unacceptable RISK at low power settings. For example, with argon beam COAGULATION, the risk of gas embolism rises if there is insufficient HF power to produce a rapid, impermeable eschar on the target tissue.
10. A warning that interference produced by the operation of HF Generator may adversely influence the operation of other electronic equipment. For PATIENTS with cardiac pacemakers or other active implants, a possible HAZARD exists because interference with the action of the active implant may occur, or the active implant may be damaged. In case of doubt, qualified advice should be obtained. Cardiologist should be consulted prior to using Sphincterotome in these patients.
11. Skin-to-skin contact should be avoided (for example between the patient's arms and body) by way of dry cloth or gauze. Monitoring electrodes should be placed as far from the surgical area as possible. Needle monitoring electrodes are not recommended. Avoid incidental contact between Active Cords and the patient's body, or any other electrodes. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the procedure.
12. Consult the neutral electrode manufacturer about the proper grounding of the patient. It is recommended that a monitoring neutral electrode be used, if a contact quality monitor is available, or built into the generator. The entire area of the neutral electrode should be attached reliably to the patient's body, and as close to the operating field as possible. The patient should not come into contact with metal parts or objects that may be grounded to earth. The use of antistatic sheeting is recommended for this purpose. Ensure that a good return path to the electrosurgical unit is maintained.
13. When HF Generator and physiological monitoring equipment are used simultaneously on the same PATIENT, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended.
In all cases, monitoring systems incorporating HF current limiting devices are recommended.

- The PATIENT leads should be positioned in such a way that contact with the PATIENT or other leads is avoided.

Apparent low output or failure of the HF Generator to function correctly at the normal operating settings may indicate faulty application of the NEUTRAL ELECTRODE or poor contact in its connections. In this case, the application of the NEUTRAL ELECTRODE and its connections should be checked before selecting a higher output power.

PREPARATION

- Remove the sphincterotome from the package.
Note: Ensure the pre-curved mandrel portion stays with the guidewire lumen of the sphincterotome.
- Remove the metal mandrel out of the sphincterotome distal tip (Be careful to keep pre-curved shape of the tip)
- Inspect the sphincterotome for proper deflection and for any damage, such as kinks.

Precaution: Kinks in the catheter will hinder injection capability. Do not use the sphincterotome if any defect is found during inspection. Please notify manufacture and return for replacement.

Precaution: Distal tip deflection greater than 90 degrees is not recommended.

WARNING: Distal tip deflection greater than 90 degrees may cause the cutting wire to break.

Precaution: Do not pre-active the sphincterotome. Pre-activating the cutting wire prior to use may cause premature cutting wire fatigue and may compromise the cutting wire's integrity.

- Ensure that the patient is properly grounded per the HF Generator manufacturer's instruction.
- If contrast media will be injected through the injection lumen, attach a syringe $\leq 20\text{ml}$ (20cc) to the proximal luer hub and inject some saline to check the injection lumen.
- The sphincterotome is now ready to use.

DIRECTIONS FOR USE

Precaution: It is recommended that the HF Generator be turned off during insertion of the sphincterotome Advance the sphincterotome using short deliberate 2-3cm movements to prevent kinking in the catheter.

WARNING: It is not allowed to open the foot switch during insertion of the sphincterotome.

- Ensure the Guidewire Locking Device is firmly seated onto the endoscope handle against the side of the working channel port; Insert the sphincterotome into the endoscope with a working channel of 3.2mm or higher.

2. Direct the cutting wire position to the desired orientation by rotating the handle either left or right.
3. Once the sphincterotome is at the desired position within the biliary system, use contrast medium to perform cholangiogram with fluoroscopy and/or x-ray or the guidewire to determine proper position within the biliary system.

WARNING: Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and endoscope while electrical current is applied. This may cause a short circuit, which may result in patient injury, broken cutting wire, and/or damage to the endoscope.

4. Prior to performing sphincterotome, make sure the sphincterotome cutting wire is in proper position when viewed endoscopically.

Precaution: If the guide wire is to be removed during sphincterotomy, turn off the power on the electrosurgical unit before removal.

5. Using the appropriate power setting recommended by the HF Generator manufacturer (30 to 70 watts), activate the cutting wire to perform sphincterotomy.
6. Upon completion of the sphincterotomy, turn off the HF Generator, disconnect the active cord, and then remove the sphincterotome.
7. To remove the sphincterotome, separate the guidewire from the open channel on the wall of the sphincterotome's catheter. Lock the guidewire on the Guidewire Locking Device, withdraw the sphincterotome until the end of open channel emerges in the device entrance of Guidewire Locking Device, unlock Guidewire and remove the Sphincterotome from endoscope completely
8. If the device in use appeared wire broken, the head severely deformed or other damage, please do not use and replace of device. Discard the sphincterotome according to local laws and/or regulations.
9. If no other device exchange is performed, unlock the guide wire from the guidewire locking device and remove from the working channel of the endoscope.

WARNINGS

Do not modify this device in any way, form or matter.

EMC CONDITIONAL

Guidance and manufacturer's declaration - electromagnetic emissions - for all ME EQUIPMENT and ME SYSTEM.

WARNING:

1. The Operator should not use the system and should inform the customer service, if the monopolar cutting is lost or degraded due to EM disturbances.

2. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the sphincterotome, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 1 Electromagnetic emission level

Electromagnetic emission	
Electromagnetic requirements of this RF generator are given below and it is the responsibility of end user to meet these requirements.	
Emission test	Compliance
CISPR 11 Conducted emission	Group 1, Class A*
CISPR 11 Radiated emission	
IEC61000-3-2 Harmonic emission	N/A
IEC61000-3-3 Voltage fluctuation / flickering emission	N/A
* The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.	

Table 2 ENCLOSURE PORT

Phenomenon	Basic EMC standard	Immunity compliant levels
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Table 4
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

Table 3 PATIENT coupling PORT

Phenomenon	Basic EMC standard	Immunity compliant levels
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

Table 4 Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	28
710	704 – 787	LTE Band 13,17	Pulse modulation 217 Hz	9
745				
780				
810	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3,4, 25; UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 – 5800	WLAN 802.11a/n	Pulse modulation 217 Hz	9
5500				
5785				

STORAGE

Store in a clean and dry environment. Do not store in direct sunlight. Do not use if labeling is incomplete or illegible. Do not use if package is damaged.

Transportation and environment	Ambient temperature	-35–60 ° C (-31–140 ° F)
	Relative humidity	0–85%
	Atmospheric pressure	700–1060 hPa (0.7–1.1 kgf/cm ²) (10.2–15.4 psia)
Operating environment	Ambient temperature	10–40 ° C (50–104 ° F)
	Relative humidity	30–85%
	Atmospheric pressure	700–1060 hPa (0.7–1.1 kgf/cm ²) (10.2–15.4 psia)

PRODUCT DISPOSAL

After use, dispose of produce and packaging in accordance with hospital, administrative and/or local government policy.

LIMITED WARRANTY AND DISCLAIMERS

1. Limited Warranty to Buyer. Micro-Tech USA warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech USA and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech USA's literature are meant to generally describe the products and do not constitute any express warranties. In the event that Micro-Tech USA gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech USA's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech USA in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

2. Disclaimer and Release. **THE WARRANTIES, OBLIGATIONS, AND LIABILITIES OF Micro-Tech USA AS SET FORTH HEREIN ARE EXCLUSIVE. BUYER HEREBY WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY BUYER, INCLUDING, BUT NOT LIMITED TO: (1) ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND (2) ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING, OR USAGE OF TRADE.**
3. Implied Warranties. The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech USA to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.
4. Limitation of Liability. **EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, Micro-Tech USA'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO: (a) THE REPLACEMENT OF THE PRODUCTS OR THE RE-SUPPLY OF EQUIVALENT PRODUCTS; (b) THE REPAIR OF THE PRODUCTS OR PAYMENT OF THE COST OF REPAIRING THE PRODUCTS; or (c) PAYMENT OF THE COST OF REPLACING THE PRODUCTS OR ACQUIRING EQUIVALENT PRODUCTS. MICRO-TECH USA SHALL HAVE NO OBLIGATION OR LIABILITY, WHETHER ARISING IN CONTRACT (INCLUDING WARRANTY), TORT (INCLUDING ACTIVE, PASSIVE, OR IMPUTED NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY) OR OTHERWISE, FOR ANY SPECIAL, CONSEQUENTIAL, PUNITIVE, INCIDENTAL, OR INDIRECT DAMAGES, OR FOR LOSS OF USE, LOSS OF REVENUE, LOSS OF BUSINESS, LOST PROFIT, OR OTHER FINANCIAL LOSS ARISING OUT OF OR IN CONNECTION WITH ANY PRODUCT OR OTHER GOODS OR SERVICES FURNISHED BY MICRO-TECH USA, EVEN IF MICRO-TECH USA WAS AWARE OF THE POSSIBILITY OF SUCH DAMAGES OR LOSS.**

SYMBOLS

	Manufacturer		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Do not reuse		Do not resterilize
	Keep away from sunlight		Keep dry
	Use by date		Date of manufacture
	Catalogue number		Batch code
	Refer to instruction manual/booklet		Do not use if package is damaged
	Sterilized using ethylene oxide		Caution
	Temperature limitation		Humidity limitation
	Type BF applied part		Contents
	Not made with natural rubber latex		Open here
	Sterile barrier system/ sterile packaging		

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech USA Inc.