# **EdgeHog™ Polypectomy Snare**

## Instructions for Use

# **Asa para polipectomía EdgeHog™**

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# Contents

---

EdgeHog™ Polypectomy Snare Instructions for use.....	1
Asa para polipectomía EdgeHog™ Instrucciones de uso.....	13

Rx only

## **IMPORTANT INFORMATION**

Caution: Federal law restricts this device to sale by or on the order of a physician. Read all instructions carefully before use. They contain essential information on using this device safely and effectively. Keep these instructions in a safe, accessible location, as you may need to refer to them again. If you have any questions or comments about any information in these instructions, please contact Micro-Tech.

## **INTENDED USE**

1. The EdgeHog™ Polypectomy Snare are used endoscopically in the removal of diminutive polyps, sessile polyps, and pedunculated polyps within the GI tract.
2. Single or multiple gastrointestinal polyps with diameter less than 2cm;
3. Non-infiltration type polyp(s) confirmed by histopathology and multiple polyps with quantity less than 30 pieces;
4. Early stage cancer only confined in mucous layer.

## **CONTRAINDICATIONS**

Contraindications include, but may not be limited to:

1. Poor physical fitness, with serious heart or/and lung disease, and can't tolerate endoscopy and endoscopic treatment;
2. Bleeding tendency, an extended bleeding and prolonged coagulation time, or thrombocytopenia or lack of prothrombin which couldn't be treated
3. A polyp / polyps with a large base more than 2cm for gastric poly and 1.5cm for normal polyp.
4. A polypoid cancer which has infiltrated into tissue or organ and has deteriorated;
5. Patients with implanted pacemaker or implanted with a metal valve / valves (It is classified as relative contraindication);

6. Diabetes, regardless of whether blood sugar is normal (poor healing capacity) (It is classified as relative contraindication);
7. Uncooperative patients or family members.

## **POTENTIAL COMPLICATIONS**

---

- Immediate or delayed hemorrhages, perforation, trans-mural burn, thermal injury, respiratory depression or apnea, arrhythmia or cardiac arrest;
- Infection, septicaemia, etc.;
- Complications which may exist and are not known or observed.

## **WARNINGS**

---


1. The product is intended for single use only! DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed .
2. Do not use this instrument for any purpose other than its intended use.
3. The product is only intended for adult populations.
4. This device is not made with natural rubber latex.
5. Patients should be informed of the potential risks and complications, which may lead to injury, illness or death of the patient.
6. The instrument is intended for use under the direct supervision of a suitably trained physician only. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
7. Confirm that the endoscopy view is clear before use. Do not insert the instrument into the endoscope unless you have a clear endoscopic

field of view. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. Damage to the endoscope and/or the instrument may also occur.

8. When applied to a patient with a pacemaker implanted, the instrument may cause malfunction or failure of the pacemaker, seriously affecting the patient. Consulting cardiologist or the manufacturer of the pacemaker is recommended before instrument application.
9. When using the instrument in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.
10. When using an electrocardiograph or other physiological monitoring equipment simultaneously with the instrument, any monitoring electrodes should be placed as far away as possible from the electrodes used with the electrosurgical unit. Needle monitoring electrodes should not be used. Physiological monitoring equipment incorporating high-frequency current limiting device is recommended.
11. When the endoscope is used with the device, the leakage currents may be relatively increased. Ensure a proper path from patient return electrode to electrosurgical unit is connected and maintained throughout the procedure.
12. Before use, please check the insertion patients body parts, ensure that no sharp edge.
13. Do not touch energized handle components (e.g. the high-frequency applied mandrel) and do not contact with conductive object as shown in Fig.1 and Fig.2 below.

**Fig.1**

**Fig.2**

14. Operation without strictly following instructions from electro-surgical unit manufacturer could cause the Snare invalid, patient or user injury.
15. The applicable electro-surgical unit should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the electro-surgical devices should be observed to verify normal operation in the configuration in which it will be used.
16. The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the electro-surgical unit manufacturers replacement parts for internal components, may result in increased emissions or decreased immunity of the electro-surgical unit.
17. Since supporting the use of Olympus endoscopic system, Olympus endoscopic system's light source may affect Polypectomy Snare's temperature to exceed 41°C, to minimize the risk of injury, please ensure the correct use of it according to Olympus endoscopic system and the manufacturer's instructions.

**【 Product Name 】** EdgeHog™ Polypectomy Snare

**【 Packaging 】** Packed in Pouch

**【 Production Date 】** See packaging

**【 Sterilization 】** Sterilized by EO (ethylene oxide) gas

**【 Shelf Life 】** See packaging

**【 Compatible Working Channel 】**  $\cong \phi$  2.8 mm

**【Applicable Electro-surgical Unit】**


High-frequency Electro-surgical unit that which is legal listing in USA is recommended such as CONMED Electro-surgery 60-8200-230 and its matched 474-L active cord.

**【Rated High-Frequency Voltage】**

DO NOT use higher repeated peak voltage than the Cut -1200Vp ( 2400Vp-p)

## STRUCTURE

The Polypectomy Snare consists of a loop assembly which can be extended and retracted from the snare's outer tube by manipulating the handle assembly. ( Refer to Fig.3 ).

**Fig.3 Structural Sketch of Sterile Snare with type II handle**

1. Oval Snare Loop 2. Outer Tube 3. Conductive Plug 4. Handle assembly

**Note:** This is the sketch map of one kind of loop, the other loop shapes are not shown here.

## PREPARATION

1. Reference the product label and choose the appropriate device.
2. Contents supplied STERILE.
3. Inspect the package before use for any damage. Do not use if package is damaged. Open the package carefully after verifying the shelf life.
4. Carefully remove the device from its packaging and uncoil it. Do NOT use excessive force as this may damage the device and affect performance.
5. Before use, please check the insertion patients body parts, ensure that no sharp edge.
6. If this device shows any signs of damage, do not use. Do not attempt to repair a nonfunctional or damaged device.
7. Before use, operate the slider back and forth to observe whether the

Snare loop extends and retracts smoothly. The Snare loop shall be free of distortion, deformation, frayed wires, disconnections and sharp protrusions. Do not use if any abnormality is found.

8. Before use, remove inflammable gas which is accumulated in human body in case of the risk of explosion.
9. Attach the negative plate of the high-frequency electrosurgical unit on the patient's vascular/muscle-rich part .
10. Inspect active cord.It must be free of kinks,bends,breaks and exposed wires to allow for accurate transfer of current.If an abnormality is noted,do not use that active cord.

**NOTE:**The active cord should be accordance with the applicable electrosurgical device accompanying documents.

---

## INSTRUCTIONS FOR USE

Hot Snare Use:

1. Retract the handle of the EdgeHog™ Polypectomy Snare and confirm that the snare loop is fully retracted into the catheter prior to inserting the catheter into the endoscope.
2. Advance the catheter through the endoscope channel until the lesion is clearly identified and targeted. Snare loop can be rotated .
3. With electrosurgical unit off, securely connect the active cord to the snare handle and electrosurgical unit.
4. Turn the electrosurgical unit ON.

NOTE: The device should be installed and put into service according to the EMC information provided in the electrosurgical unit accompanying documents.

NOTE: In strict accordance with the applicable electrosurgical device operational mode and Input feature during operation.

Take the CONMED Electrosurgery 60-8200-230 as an example:

Operational mode: intermittence 15s open/30s close;

The rated Input:100-240V,4A max.50/60Hz.

5. Provide a generous submucosal lift to fully separate the lesion from the muscularis propria.

NOTE: Do not use in Underwater EMR (UEMR) without submucosal lift.

6. Operate the handle to extend the snare loop and indent snare into the target tissue. Maintain visualization of the distal part of the snare loop. Close snare completely on the tissue and ensure desired tissue is captured. If necessary, regrip on the tissue to entrap less. Press the foot switch to active output. Finish pulling the slider and resect the target tissue. Do not apply excessive force when snaring the tissue. This could cause patient injury, such as punctures, hemorrhages or mucous membrane damage or thermal injury.
7. Retract the snare prior to removing the instrument from the endoscope
8. Upon completion of the procedure, disconnect the active cord from the snare handle and dispose of the device according to the institutional guidelines for biohazard medical waste and local regulatory requirements.

## EMC CONDITIONAL

Guidance and manufacturer's declaration - electromagnetic emissions - for all ME EQUIPMENT and ME SYSTEM.

Table 1: Guidance and manufacturer's declaration – electromagnetic emissions		
The Polypectomy Snare is intended for use in the electromagnetic environment specified below. The customer or the user of the Polypectomy Snare should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Polypectomy Snare uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class A	The Polypectomy Snare is suitable for use in all establishments other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000 -3-2	A	
Voltage fluctuations/ Flicker emissions IEC 61000 -3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity - for all ME EQUIPMENT and ME SYSTEM.

Table 2:Guidance and manufacturer's declaration – electromagnetic immunity			
The Polypectomy Snare is intended for use in the electromagnetic environment specified below. The customer or the user of the Polypectomy Snare should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV lines to lines ±2 kV lines to earth	±1 kV lines to lines ±2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) For 0,5 cycle 40% UT (60% dip in UT) For 5 cycle 70% UT (30% dip in UT) For 25 cycle <5% UT (>95% dip in UT) For 5 cycle	<5% UT (>95% dip in UT) For 0,5 cycle 40% UT (60% dip in UT) For 5 cycle 70% UT (30% dip in UT) For 25 cycle <5% UT (>95% dip in UT) For 5 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of Polypectomy Snare requires continued operation during power mains interruptions, it is recommended that the Polypectomy Snare be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
Note: UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity for ME EQUIPMENT and ME SYSTEM that are not LIFE-SUPPORTING.

Table 3: Guidance and manufacturer's declaration – electromagnetic immunity			
The Polypectomy Snare is intended for use in the electromagnetic environment specified below. The customer or the user of the Polypectomy Snare should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6			Portable and mobile RF communications equipment should be used no closer to any part of the SYSTEM, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80MHz	3V	$d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80MHz to 800MHz  $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800MHz to 2.5GHz
	3 V/m 80MHz to 2.5GHz	3V/m	Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.  Interference may occur in the vicinity of equipment marked with the following symbol: 

Recommended separation distances between portable and mobile  
RF communications equipment and the ME EQUIPMENT and ME  
SYSTEM – for ME EQUIPMENT and ME SYSTEM that are not  
LIFE-SUPPORTING.

Recommended separation distances between portable and mobile RF communications equipment and the Polypectomy Snare			
The Polypectomy Snare is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Polypectomy Snare can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters ) and the Polypectomy Snare as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80 MHz $d = [\frac{3.5}{V_1}] \sqrt{P}$	80MHz to 800MHz $d = [\frac{3.5}{E_1}] \sqrt{P}$	800MHz to 2.5GHz $d = [\frac{7}{E_1}] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

## **STORAGE**

The product should be stored in a cool, dry, clean, well-ventilated, non-corrosive gas environment.

Do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

## **PRODUCT DISPOSAL**

After use, this product may be a potential biohazard. Handle and dispose of in accordance with hospital, local and administrative laws and regulations.

## **Limited Warranty and Disclaimers :**

---

1. Limited Warranty to Buyer. Micro-Tech USA warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech USA and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech USA's literature are meant to generally describe the products and do not constitute any express warranties. In the event that Micro-Tech USA gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech USA's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech USA in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.
2. Disclaimer and Release. THE WARRANTIES, OBLIGATIONS, AND LIABILITIES OF Micro-Tech USA AS SET FORTH HEREIN ARE EXCLUSIVE. BUYER HEREBY WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY BUYER, INCLUDING, BUT NOT LIMITED TO: (1) ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND (2) ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING, OR USAGE OF TRADE.
3. Implied Warranties. The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech USA to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.
4. Limitation of Liability. EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, Micro-Tech USA'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO: (a) THE REPLACEMENT OF THE PRODUCTS OR THE RE-SUPPLY OF EQUIVALENT PRODUCTS; (b) THE REPAIR OF THE PRODUCTS OR PAYMENT OF THE COST OF REPAIRING THE PRODUCTS; or (c) PAYMENT

OF THE COST OF REPLACING THE PRODUCTS OR ACQUIRING EQUIVALENT PRODUCTS. MICRO-TECH USA SHALL HAVE NO OBLIGATION OR LIABILITY, WHETHER ARISING IN CONTRACT (INCLUDING WARRANTY), TORT (INCLUDING ACTIVE, PASSIVE, OR IMPUTED NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY) OR OTHERWISE, FOR ANY SPECIAL, CONSEQUENTIAL, PUNITIVE, INCIDENTAL, OR INDIRECT DAMAGES, OR FOR LOSS OF USE, LOSS OF REVENUE, LOSS OF BUSINESS, LOST PROFIT, OR OTHER FINANCIAL LOSS ARISING OUT OF OR IN CONNECTION WITH ANY PRODUCT OR OTHER GOODS OR SERVICES FURNISHED BY MICRO-TECH USA, EVEN IF MICRO-TECH USA WAS AWARE OF THE POSSIBILITY OF SUCH DAMAGES OR LOSS.

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**Micro-Tech (Nanjing) Co., Ltd.**
**Micro-Tech USA Inc.**