

DEPENDABLE & PRECISE ENDOSCOPIC ULTRASOUND NEEDLES

Whether requiring an EUS solution for diagnostic and therapeutic intervention or for histological sampling, Micro-Tech Endoscopy provides high quality, affordable clinical options to help with consistent sampling and accurate diagnosis from EUS tissue acquisition.

The Trident™ FNB Needle with its unique tri-tip core needle design achieves maximized yield and enables consistent histological tissue acquisition.

The Areus™ FNA Needle features an elongated dual edged bevel tip to aid in precise penetration into the target area delivering a consistent yield for cytology.



NEVER
COMPROMISE
ON QUALITY

DRAMATICALLY
IMPROVE YOUR
BOTTOM LINE

THROW AWAY
CONTRACTS
FOREVER

RELIABLE
SUPPLY
PARTNERSHIP

KEY BENEFITS

COBALT-CHROMIUM (Co-Cr) NEEDLE CONSTRUCTION

The Trident and Areus endoscopic ultrasound needles are made from the Co-Cr alloy. Co-Cr maintains its shape and sharpness while providing durability, flexibility and reliable target accessibility.

ECHOGENICITY

It is always important to maintain tip visibility for accurate guidance within the target area. The surface of the Trident and Areus needles feature a laser etch V pattern, that provides good visibility under ultrasound.

SIMPLE PRECISE HANDLE CONTROLS

The Trident and Areus needles offer an ergonomic handle design with a smooth slide action during tissue sampling. Push button adjusters for the needle and sheath length provide comfortable and precise one-handed adjustment and locking control.

SPECIFICATIONS

TRIDENT® ENDOSCOPIC ULTRASOUND FINE NEEDLE BIOPSY

Order Number	Henry Schein Item Number	Description	Needle Size	Sheath Diameter (mm)	Minimum Channel Size (mm)	Package Units
ES34611	TBD	Trident FNB Needle	25ga	1.8	2.8	Each
ES34591	TBD	Trident FNB Needle	22ga	1.8	2.8	Each

AREUS® ENDOSCOPIC ULTRASOUND ASPIRATION NEEDLE

Order Number	Henry Schein Item Number	Description	Needle Size	Sheath Diameter (mm)	Minimum Channel Size (mm)	Package Units
ES34621	TBD	Areus FNA Needle	25ga	1.8	2.8	Each
ES34601	TBD	Areus FNA Needle	22ga	1.8	2.8	Each

- Working Length: 137.5cm to 141.5cm, adjustable
- Needle Length: 0 to 8cm, adjustable
- Includes Vacuum Lock 20cc syringe and one-way stopcock

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



EXCEPTIONAL
ENDOSCOPIC
SOLUTIONS.

2855 BOARDWALK DRIVE
ANN ARBOR, MI 48104
877.552.4027

MTENDOSCOPY.COM

Customercare@mtendoscopy.com
EUS_MPW43000 Rev.3