



ClearPetra™ Scope/Sheath Cross Reference List

Please use the following cross reference list to determine the minimum size ClearPetra™ sheaths you need for the scopes you use. **Note: There must be a 3Fr difference between your scope outer diameter (OD) and the ClearPetra™ sheath inner diameter (ID).*

Flexible Ureteroscopy Instruments

Brand	Model	OD (Fr)	Minimum ID of ClearPetra
Richard Wolf	BOA	8.7Fr	12Fr CP sheath
Richard Wolf	Cobra	9.9Fr	13Fr CP sheath
Olympus	URF-V	9.9Fr	13Fr CP sheath
Olympus	URF-P6	7.95Fr	11Fr CP sheath
Olympus	URF-V2	8.4Fr	12Fr CP sheath
Olympus	URF-P5	8.4 Fr	12Fr CP sheath
Olympus	URF-P3	8.4 Fr	12Fr CP sheath
Olympus	URF-P2	9.9Fr	13Fr CP sheath
Gyrus - ACMI	MR-6A	6.9Fr	10Fr CP sheath
Karl Storz	Flex-Xc	8.5Fr	12Fr CP sheath
Karl Storz	Flex-X2S	7.5Fr	11Fr CP sheath
Stryker	U-500	6.9Fr	10Fr CP sheath
Boston	LithoVue	9.5Fr	13Fr CP sheath
Pusen		9.5Fr	13Fr CP sheath

Semi-rigid Ureteroscopy Instruments

Brand	Model	OD (Fr)	Minimum ID of ClearPetra
Wolf	8701.543	4.5Fr-6.5Fr	10Fr CP sheath
Wolf	8702.5X4	6Fr-7.5Fr	11Fr CP sheath
Wolf	8703.5X4	8Fr-9.8Fr	13Fr CP sheath
Wolf	8708.518	6.5Fr-8.5Fr	12Fr CP sheath
Wolf	8708.534	6.5Fr-8.5Fr	12Fr CP sheath
Storz	27000K/L	6.5Fr-7Fr-9.9Fr	13Fr CP sheath
Storz	27010K/L	7Fr-8.4Fr-9.9Fr	13Fr CP sheath
Storz	27002KP	7.3Fr-8Fr	11Fr CP sheath
Olympus	WA29040/1A	6.4Fr-7.8Fr	12Fr CP sheath
Olympus	WA29042A	8.6Fr-9.8Fr	13Fr CP sheath
Olympus	WA29048/9A	6.4Fr-7.8Fr	12Fr CP sheath
Olympus	WA02943/4/6A	7.5Fr	11Fr CP sheath
Stryker	505-880-3/430	6.9Fr-10Fr	13Fr CP sheath
Sholly	41.0611a.wol	7.5Fr-9.8Fr	13Fr CP sheath
Sholly	41.0609a.wol	7.5Fr-9.8Fr	13Fr CP sheath
Sholly	41.0615a.wol	6.0Fr-7.5Fr	12Fr CP sheath
Sholly	41.0614a.wol	6.0Fr-7.5Fr	12Fr CP sheath



ClearPetra™ scope-sheath cross reference list (continued)

Nephrostomy list

Brand	Model	OD (Fr) and length	Minimum ID of ClearPetra
Storz	MIP M,27830 kak	12Fr, 220mm	16Fr in 13/15cm lengths
Storz	MIP S/XS 27820 kak	7.5Fr, 240mm	11Fr in 13/15/17cm lengths
Wolf	Lahme,8968421	12Fr, 225mm	16Fr in 13/15cm lengths
Olympus	Mini	11Fr, 220mm	14Fr in 13/15cm lengths
Scholly	Nephroscope 41.0634a	13.5Fr, 220mm	18Fr in 13/15cm lengths
Scholly	Mini 41.0635a	13.5Fr, 320mm	18Fr in 13/15/17/21cm lengths

Please note: This list is subject to change and may not list all available scopes.
Please be sure to confirm the measurements of your scope when determining sheath size to ensure you will have a 3Fr difference between the scope and sheath.