

# EFFECTIVE LIGATION MULTIPLE BAND LIGATORS



Endoscopic banding has been demonstrated as an effective and simple treatment for internal hemorrhoids and varices. Hemorrhoid banding is readily performed in an ambulatory setting and can result in less pain and have a shorter recovery period than a surgical hemorrhoidectomy.

The Multi-Band Ligator Sets from Micro-Tech offer cost-effective solutions for the treatment of hemorrhoids and varices during routine endoscopy procedures.

## KEY BENEFITS

### EXCELLENT VISIBILITY

A clear ligating barrel with rear mounted bands optimizes the endoscopic view while maintaining suction.

### BAND DESIGN

The square band design provides consistent shape and reliable retention. The hemorrhoid set comes with four pre-loaded latex-free bands for routine colonoscopies. The upper ligator comes pre-loaded with six latex-free bands and can be used in upper and lower endoscopic procedures.

### USER-FRIENDLY SYSTEM

The handle unit tie hole enables quick and easy loading of the pull wire. The alert band notifies the operator when only one band remains.

## SPECIFICATIONS

Order Number	Henry Schein Item Number	Description	Minimum Working Channel (mm)	Working Length (cm)	Package Units
HL32161	131-9042	Endoscopic Hemorrhoid Multi-Band Ligator Set, 4 Latex Free Bands, compatible with 9.4-13 mm scope diameter	2.8	220	1/Box
ML32651	138-5348	Endoscopic Multi-Band Ligator Set, 6 Latex Free Bands, compatible with 9.4-13 mm scope diameter	2.8	145	1/Box