

PERSONAL PROTECTION DISPOSABLE MEDICAL FACE MASK



COVID-19 has brought unique challenges to our world today and Micro-Tech Endoscopy is committed to supporting our health care professionals with new products, valuable education opportunities and access to much-needed Personal Protection Equipment (PPE).

The 3-ply Medical Face Mask is intended for use in healthcare settings to help prevent the spread of infection or illness during the COVID-19 pandemic. The product has not been FDA cleared or approved and is authorized by FDA under an EUA for use as source control by the general public as well as by HCP in healthcare settings.

- Ear loop style, 175mm x 95mm single use
- Used for source control in medical settings
- Normal donning and doffing procedures apply
- Not to be used for particulate filtration
- Do not use in high-risk aerosol generating procedures

DISPOSABLE MEDICAL FACE MASK

Order Number	Description	Package Units
NKY-KZ-175X95	Medical Face Mask	50/Box