

THE JACKET™ DISPOSABLE BIOPSY VALVE



The Jacket™ Disposable Biopsy Valve is used to cover the opening to the biopsy/suction channel of gastrointestinal endoscopes.

Its unique design offers improved valve performance by minimizing leakage and helping to preserve insufflation to maintain adequate visualization.

- Double valve system to minimize leakage and prevent spraying
- Improved seal with three points of contact
- Internal funnel seal hugs and flexes with inserted accessories without gaping
- Compatible with Olympus, Pentax and Fuji GI endoscopes



SPECIFICATIONS

THE JACKET™ DISPOSABLE BIOPSY VALVE

Order Number	Henry Schein Item Number	Description	Package Units
QL70021	138-6532	Jacket Disposable Biopsy Valve, Olympus/Fujinon, Non-Sterile (Blue)	Box of 100
QL70051	138-6534	Jacket Disposable Biopsy Valve, Pentax Medical, Non-Sterile (Clear)	Box of 100

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.