

## ABC<sup>™</sup> 3-STEP INFLATION BALLOONS

For a simple, affordable solution to balloon dilation, choose the ABC Dilation Balloon. Manufactured from a high-quality non-compliant nylon material, the ABC Balloon offers the radial strength required to accurately and consistently dilate strictures.

The balloon can be inflated in 3 steps, eliminating the need for multiple balloons in the same procedure. It offers an economical solution versus typical radial expansion balloons.

The balloon is available as a flexible fixed wire esophageal 8 cm balloon or a fixed wire esophageal, pyloric and colonic 5.5 cm balloon, both featuring rapid inflation/deflation performance.



## KEY FEATURES

### CHOICE OF MATERIAL

The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilatation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

### ROUNDED PROXIMAL SHOULDER

The combination of a rounded proximal shoulder and high-quality transparent nylon material facilitates visualization of the stricture during dilation.

### BALLOON LENGTH

Choose from a balloon length of 5.5 cm, which can be used as a universal device to dilate esophageal, pyloric or colonic strictures, or an 8 cm balloon depending on clinical need.

## SPECIFICATIONS

### ABC DILATION BALLOON

Order Number	Henry Schein® Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24165	139-6667	ABC Dilation Balloon	6	5.5	230	2.8	2/Box
FW24175	139-6668	ABC Dilation Balloon	8	5.5	230	2.8	2/Box
FW24185	139-6669	ABC Dilation Balloon	10	5.5	230	2.8	2/Box
FW24195	139-6670	ABC Dilation Balloon	12	5.5	230	2.8	2/Box
FW24205	139-6671	ABC Dilation Balloon	15	5.5	230	2.8	2/Box
FW24215	139-6672	ABC Dilation Balloon	18	5.5	230	2.8	2/Box
FW24225	139-6673	ABC Dilation Balloon	20	5.5	230	2.8	2/Box

### ABC ESOPHAGEAL DILATION BALLOON

Order Number	Henry Schein® Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	128-5660	ABC Esophageal Dilation Balloon	6	8	180	2.8	2/Box
FW24085	128-5661	ABC Esophageal Dilation Balloon	8	8	180	2.8	2/Box
FW24095	128-5662	ABC Esophageal Dilation Balloon	10	8	180	2.8	2/Box
FW24105	128-5663	ABC Esophageal Dilation Balloon	12	8	180	2.8	2/Box
FW24125	128-5664	ABC Esophageal Dilation Balloon	15	8	180	2.8	2/Box
FW24145	128-5665	ABC Esophageal Dilation Balloon	18	8	180	2.8	2/Box
FW24155	128-5666	ABC Esophageal Dilation Balloon	20	8	180	2.8	2/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
[customer@mtendoscopy.com](mailto:customer@mtendoscopy.com)

MPW24000 Rev.4