

CLEARLY, THE ULTIMATE POLYP TRAP PROTRAP LUMA



The ProTrap Luma™ in-line polyp trap and specimen tray consist of a collection hose, a trap body, a tray, and a glow stick. Its unique design offers the potential for touch-less sample retrieval. The specimen tray can be removed and placed directly into the specimen jar, minimizing handling and the associated risks. Additional trays are available as needed, and the trap is complemented with a light source for ease of use.

- Bright and long-lasting glow stick helps aid in retrieval visualization
- The entire tray fits inside the specimen jar
- Built-in magnifying window for easy visualization of specimen



SPECIFICATIONS

PROTRAP LUMA [Ⓢ]

Order Number	Henry Schein Item Number	Description	Package Units
TP54411	TBD	ProTrap Luma In-Line Polyp Trap	25/Box
TP54421	TBD	ProTrap Luma In-Line Polyp Trap Specimen Tray	50/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.