

# CLEAR DIFFERENTIATION BLUEBOOST SUBMUCOSAL INJECTABLE SOLUTION

BlueBoost™ is a lifting agent that lasts. Once injected into the submucosal space, BlueBoost adds efficiency to GI procedures by separating the layers of the intestine and helping to delineate the margin around the lesion. This premixed deep sky-blue viscous solution is designed to assist with removal of polyps, adenomas and other lesions.

- Comes in prefilled syringes that are ready to use
- Outlasts saline and creates a reliable lift that will save you time



## SPECIFICATIONS

**BLUEBOOST** STERILE 

Order Number	Description	Package Units
SL62521	BlueBoost™ Submucosal Injectable Solution, 10ml	12/Box

## CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.