


ESD SPECIMEN TRAY PROTECT TISSUE SAMPLES



This convenient specimen container has everything needed to mount tissue specimens before sending them to pathology. The tabs on the lid help seal the box so tissue is protected while it is transported. A tube of push pins is included in the box for mounting the tissue to the specimen pad. Markings on the blue pad help track orientation and include a measurement grid for sizing.



SPECIFICATIONS

ESD SPECIMEN TRAY ⁽²⁾

Order Number	Description	Package Units	Dimensions
QB55011	ESD Specimen Tray with pins	6/Box	5.5" x 5.5" x 2.75"

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.