

## IDEAL BITE TWIN-WIRE BIOPSY FORCEPS

Twin-wire mechanism and precision stamping technology for assured tissue sampling.

- Pivoting jaws to aid in sampling where space may be limited
- Elongated jaw for greater tissue retention
- With-needle option for multiple bites in a single pass
- Entry - Exit indicator markings on sheath
- Flexible and pushable smooth catheter



# KEY FEATURES

## TWIN-WIRE DESIGN

The twin-wire design allows the forcep head to rotate radially, while providing suitable closure force to the jaws, which may improve sampling ability and complete jaw closure.

## ERGONOMIC HANDLE

Three-ring comfort with cushioned lining and short throw provides improved grip and accommodates smaller hands.

## SUITABLE SAMPLE SIZE

The high yield family of forceps from Micro-Tech help the physician to achieve a suitable amount of tissue for accurate and thorough diagnosis. The large capacity, jumbo and jumbo max forceps are designed to provide optimal bite along with spike options for multiple bites in a single pass.



# SPECIFICATIONS STERILE

## TWIN-WIRE BIOPSY FORCEPS

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40286	133-4667	Large Capacity Forceps	Alligator	Blue	230	2.8	50/box
BF40306	133-4666	Large Capacity Forceps with Spike	Alligator	Blue	230	2.8	50/box
BF40461	TBD	Jumbo Biopsy Forceps	Serrated	Blue	230	2.8	10/box
BF40451	TBD	Jumbo Biopsy Forceps with Spike	Serrated	Blue	230	3.2	10/box
BF40501	TBD	Jumbo Max Biopsy Forceps	Serrated	Blue	230	3.7	10/box
BF40491	TBD	Jumbo Max Biopsy Forceps with Spike	Serrated	Blue	230	3.7	10/box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
[customer@mtendoscopy.com](mailto:customer@mtendoscopy.com)