

BREATHE-EZ BITE BLOCK PREMIUM MATERIAL AND ENHANCED DESIGN



Micro-Tech Bite Blocks offer comfort and protection during upper endoscopic procedures. The new Breathe-EZ™ Bite Block was designed with a special, softer bite plate that helps to keep the lower jaw forward and the airway open.

- Unique design that may assist in improving oxygenation during endoscopic procedures
- Available in rigid hook and loop strap and adjustable elastic strap
- Softer overlay that helps protect patient's teeth



SPECIFICATIONS

BREATHE-EZ BITE BLOCK

Order Number	Henry Schein Item Number	Description	Package Units
QL26021	TBD	Breathe-EZ Bite Block, Adult, 60Fr, with hook and loop strap	50/Box
QL26031	TBD	Breathe-EZ Bite Block, Adult, 60Fr, with elastic strap	50/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.