

## ABC 3-STEP INFLATION BALLOONS

For a simple, affordable solution to balloon dilation, choose the ABC™ Dilation Balloon. Manufactured from a high-quality non-compliant nylon material, the ABC Balloon offers the radial strength required to accurately and consistently dilate strictures.

The balloon can be inflated in 3 steps to provide gradual dilation of the stricture with a single balloon. It offers an economical solution versus typical radial expansion balloons.

The balloon is available as a flexible fixed wire esophageal 8 cm balloon or a fixed wire esophageal, pyloric and colonic 5.5 cm balloon, both featuring rapid inflation/deflation performance.



### KEY FEATURES

#### CHOICE OF MATERIAL

The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilatation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

#### ROUNDED PROXIMAL SHOULDER

The combination of a rounded proximal shoulder and high-quality transparent nylon material facilitates visualization of the stricture during dilation.

#### BALLOON LENGTH

Choose from a balloon length of 5.5 cm, which can be used as a universal device to dilate esophageal, pyloric or colonic strictures, or an 8 cm balloon depending on clinical need.

# SPECIFICATIONS

## ABC ESOPHAGEAL, PYLORIC, COLONIC DILATION BALLOON STERILE

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24165	139-6667	ABC Dilation Balloon	6	5.5	230	2.8	2/Box
FW24175	139-6668	ABC Dilation Balloon	8	5.5	230	2.8	2/Box
FW24185	139-6669	ABC Dilation Balloon	10	5.5	230	2.8	2/Box
FW24195	139-6670	ABC Dilation Balloon	12	5.5	230	2.8	2/Box
FW24205	139-6671	ABC Dilation Balloon	15	5.5	230	2.8	2/Box
FW24215	139-6672	ABC Dilation Balloon	18	5.5	230	2.8	2/Box
FW24225	139-6673	ABC Dilation Balloon	20	5.5	230	2.8	2/Box

## ABC ESOPHAGEAL DILATION BALLOON STERILE

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	128-5660	ABC Esophageal Dilation Balloon	6	8	180	2.8	2/Box
FW24085	128-5661	ABC Esophageal Dilation Balloon	8	8	180	2.8	2/Box
FW24095	128-5662	ABC Esophageal Dilation Balloon	10	8	180	2.8	2/Box
FW24105	128-5663	ABC Esophageal Dilation Balloon	12	8	180	2.8	2/Box
FW24125	128-5664	ABC Esophageal Dilation Balloon	15	8	180	2.8	2/Box
FW24145	128-5665	ABC Esophageal Dilation Balloon	18	8	180	2.8	2/Box
FW24155	128-5666	ABC Esophageal Dilation Balloon	20	8	180	2.8	2/Box

## INFLATION DEVICE STERILE

Order Number	Henry Schein Item Number	Description	Package Units
DB25011	131-4493	Micro-Tech Balloon Inflation Device 60 ml with Lock	1/Box
DB25012	TBD	Micro-Tech Balloon Inflation Device 60 ml with Lock	5/Case

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
customer@mtendoscopy.com

MPW24000 Rev. 5