


TRUE BLUE BLUE BEACON GI MARKER



- Blue Beacon™ is useful in GI procedures that require a contrast stain to temporarily mark or highlight tissue
- Can be diluted and introduced to the irrigation, sprayed or injected around the area of interest
- The 0.4% Indigo Carmine solution is supplied sterile in a convenient luer lock 5ml syringe, 5 syringes per box
- The 0.8% Indigo Carmine solution is supplied sterile in a 10 ml single use vial, 10 vials per box



SPECIFICATIONS

BLUE BEACON GI MARKER STERILE

Order Number	Henry Schein Item Number	Description	Package Units
IC62021	128-5667	Blue Beacon GI Marker 5 ml Syringe Indigo Carmine 0.4%	5/Box
IC62061	137-5185	Blue Beacon GI Marker 10 ml Vial Indigo Carmine 0.8%	10/Box

This solution is NOT FOR INTRAVENOUS OR INTRAMUSCULAR INJECTION.
It is a single use device indicated for temporary marking of the GI tract.

SPRAY CATHETER STERILE

Order Number	Henry Schein Item Number	Description	Working Length (cm)	Minimum Channel Size (mm)	Package Units
SC62241	133-5544	Spray Catheter	220	2.2	20/Box
SC62301	133-5545	Spray Catheter	250	2.8	20/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.