

NECROSECTOMY

ENDOROTOR POWERED ENDOSCOPIC DEBRIDEMENT




The EndoRotor[®] PED[™] Catheter is the first clinically validated solution for endoscopic debridement of walled-off necrosis following EUS guided drainage.

Non-thermal mechanical debridement and removal of necrotic tissue in symptomatic WOPN/WON.

Debridement of WOPN/WON under direct visualization without the need to pass the scope in and out of the collection.

360° catheter rotation to facilitate precise placement within the collection.



ENDOROTOR[®] PURGE KIT



POWERED ENDOSCOPIC
DEBRIDEMENT CATHETER

SPECIFICATIONS

ENDOROTOR POWERED ENDOSCOPIC DEBRIDEMENT (PED)

Order Number	Description	Package Units
ESS-CONSOLE	EndoRotor System	1/Box
3.2-PED-EGD-OP	Powered Endoscopic Debridement Cath, Gastroscope, 3.2 mm (Olympus/Pen)	1/Box
3.2-PED-EGD-F	Powered Endoscopic Debridement Cath, Colonoscope, 3.2 mm (Fuji)	1/Box
PED-PURGE	EndoRotor Purge Kit	10/Box
ESS-VAC-BAG	Suction Bag 1.5l	40/Box