


# STRENGTH AND CONTROL

## MULTI-STAGE DILATION BALLOONS




Choose strength and control when dilating strictures with a wire-guided or fixed wire multi-stage dilatation balloon from Micro-Tech. Multi-Stage balloons feature a premium clear Pebax<sup>®</sup> balloon with rounded shoulders for balloon endoscopy.

Multi-stage balloons allow for inflation to 3 distinct sizes to allow multi-step dilation in one balloon, providing controlled dilation and eliminating the need for multiple balloons in a procedure.

Fixed wire esophageal balloons offer fast, reliable performance.

Pre-loaded .035 wire-guided balloons assist placement in tortuous or altered anatomy.



### KEY BENEFITS

#### CHOICE OF MATERIAL

The high quality Pebax material produces a balloon with the features and performance your providers demand. Accurate three-size radial expansion, rounded shoulders, clear visibility and strength with control.

#### VISIBILITY UNDER FLUOROSCOPY

All wire-guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

#### TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire-guided balloon is provided with a pre-loaded .035 guidewire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

# SPECIFICATIONS

## WIRE-GUIDED ESOPHAGEAL, PYLORIC, COLONIC MULTI-STAGE BALLOONS STERILE

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
MB23031	125-2734	Reliant Multi-Stage Dilatation Balloon	6-7-8	5.5	230	2.8	1/Box
MB23071	125-2736	Reliant Multi-Stage Dilatation Balloon	8-9-10	5.5	230	2.8	1/Box
MB23111	125-2738	Reliant Multi-Stage Dilatation Balloon	10-11-12	5.5	230	2.8	1/Box
MB23151	125-2740	Reliant Multi-Stage Dilatation Balloon	12-13.5-15	5.5	230	2.8	1/Box
MB23191	125-2742	Reliant Multi-Stage Dilatation Balloon	15-16.5-18	5.5	230	2.8	1/Box
MB23221	125-2744	Reliant Multi-Stage Dilatation Balloon	18-19-20	5.5	230	2.8	1/Box

## FIXED WIRE ESOPHAGEAL MULTI-STAGE BALLOONS STERILE

Order Number	Henry Schein Item Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24015	129-0163	Rapide Esophageal Dilatation Balloon	6-7-8	8	180	2.8	2/Box
FW24025	129-0164	Rapide Esophageal Dilatation Balloon	8-9-10	8	180	2.8	2/Box
FW24035	129-0165	Rapide Esophageal Dilatation Balloon	10-11-12	8	180	2.8	2/Box
FW24045	129-0166	Rapide Esophageal Dilatation Balloon	12-13.5-15	8	180	2.8	2/Box
FW24055	129-0167	Rapide Esophageal Dilatation Balloon	15-16.5-18	8	180	2.8	2/Box
FW24065	129-0168	Rapide Esophageal Dilatation Balloon	18-19-20	8	180	2.8	2/Box

## INFLATION DEVICE STERILE

Order Number	Henry Schein Item Number	Description	Package Units
DB25011	131-4493	Micro-Tech Balloon Inflation Device 60 ml with Lock	1/Box
DB25012	TBD	Micro-Tech Balloon Inflation Device 60 ml with Lock	5/Case

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
[customercare@mtendoscopy.com](mailto:customercare@mtendoscopy.com)

MPW23000 Rev.5