


ABC 3-STEP INFLATION BALLOONS

For a simple, affordable solution to balloon dilation, choose the ABC™ Dilation Balloon. Manufactured from a high-quality non-compliant nylon material, the ABC Balloon offers the radial strength required to accurately and consistently dilate strictures.

The balloon can be inflated in 3 steps to provide gradual dilation of the stricture with a single balloon. It offers an economical solution versus typical radial expansion balloons.

The balloon is available as a flexible fixed wire esophageal 8 cm balloon or a fixed wire esophageal, pyloric and colonic 5.5 cm balloon, both featuring rapid inflation/deflation performance.



KEY FEATURES

CHOICE OF MATERIAL

The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilatation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

ROUNDED PROXIMAL SHOULDER

The combination of a rounded proximal shoulder and high-quality transparent nylon material facilitates visualization of the stricture during dilation.

BALLOON LENGTH

Choose from a balloon length of 5.5 cm, which can be used as a universal device to dilate esophageal, pyloric or colonic strictures, or an 8 cm balloon depending on clinical need.

SPECIFICATIONS

ABC ESOPHAGEAL, PYLORIC, COLONIC DILATION BALLOON STERILE

Order Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24165	ABC Dilation Balloon	6	5.5	230	2.8	2/Box
FW24175	ABC Dilation Balloon	8	5.5	230	2.8	2/Box
FW24185	ABC Dilation Balloon	10	5.5	230	2.8	2/Box
FW24195	ABC Dilation Balloon	12	5.5	230	2.8	2/Box
FW24205	ABC Dilation Balloon	15	5.5	230	2.8	2/Box
FW24215	ABC Dilation Balloon	18	5.5	230	2.8	2/Box
FW24225	ABC Dilation Balloon	20	5.5	230	2.8	2/Box

ABC ESOPHAGEAL DILATION BALLOON STERILE

Order Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	ABC Esophageal Dilation Balloon	6	8	180	2.8	2/Box
FW24085	ABC Esophageal Dilation Balloon	8	8	180	2.8	2/Box
FW24095	ABC Esophageal Dilation Balloon	10	8	180	2.8	2/Box
FW24105	ABC Esophageal Dilation Balloon	12	8	180	2.8	2/Box
FW24125	ABC Esophageal Dilation Balloon	15	8	180	2.8	2/Box
FW24145	ABC Esophageal Dilation Balloon	18	8	180	2.8	2/Box
FW24155	ABC Esophageal Dilation Balloon	20	8	180	2.8	2/Box

INFLATION DEVICE STERILE

Order Number	Description	Package Units
DB25011	Micro-Tech Balloon Inflation Device 60 ml with Lock	1/Box
DB25012	Micro-Tech Balloon Inflation Device 60 ml with Lock	5/Case

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTEndoscopy.com
customer@mtendoscopy.com

MPW24000 Rev. 6