

STRENGTH AND CONTROL

MULTI-STAGE DILATION BALLOONS

EXPANDS TO
3 DISTINCT SIZES

RADIOPAQUE BANDS

ROUNDED SHOULDERS

LUMINOUS LABEL

.035 COMPATIBLE

Choose strength and control when dilating strictures with a wire-guided or fixed wire multi-stage dilatation balloon from Micro-Tech. Multi-Stage balloons feature a premium clear Pebax[®] balloon with rounded shoulders for balloon endoscopy.

Multi-stage balloons allow for inflation to 3 distinct sizes to allow multi-step dilation in one balloon, providing controlled dilation and eliminating the need for multiple balloons in a procedure.

Fixed wire esophageal balloons offer fast, reliable performance.

Pre-loaded .035 wire-guided balloons assist placement in tortuous or altered anatomy.

KEY BENEFITS

CHOICE OF MATERIAL

The high quality Pebax material produces a balloon with the features and performance your providers demand. Accurate three-size radial expansion, rounded shoulders, clear visibility and strength with control.

VISIBILITY UNDER FLUOROSCOPY

All wire-guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire-guided balloon is provided with a pre-loaded .035 guidewire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

EXCEPTIONAL
ENDOSCOPIC
SOLUTIONS

Multistage Radial Expansion	
Inflation Pressure	ATM
304	3
456	4.5
608	6
Balloon O.D.	
MM	Ft
18	54
19	57
20	60

SPECIFICATIONS

WIRE-GUIDED ESOPHAGEAL, PYLORIC, COLONIC MULTI-STAGE BALLOONS STERILE

Order Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
MB23031	Reliant Multi-Stage Dilatation Balloon	6-7-8	5.5	230	2.8	1/Box
MB23071	Reliant Multi-Stage Dilatation Balloon	8-9-10	5.5	230	2.8	1/Box
MB23111	Reliant Multi-Stage Dilatation Balloon	10-11-12	5.5	230	2.8	1/Box
MB23151	Reliant Multi-Stage Dilatation Balloon	12-13.5-15	5.5	230	2.8	1/Box
MB23191	Reliant Multi-Stage Dilatation Balloon	15-16.5-18	5.5	230	2.8	1/Box
MB23221	Reliant Multi-Stage Dilatation Balloon	18-19-20	5.5	230	2.8	1/Box

FIXED WIRE ESOPHAGEAL MULTI-STAGE BALLOONS STERILE

Order Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24015	Rapide Esophageal Dilatation Balloon	6-7-8	8	180	2.8	2/Box
FW24025	Rapide Esophageal Dilatation Balloon	8-9-10	8	180	2.8	2/Box
FW24035	Rapide Esophageal Dilatation Balloon	10-11-12	8	180	2.8	2/Box
FW24045	Rapide Esophageal Dilatation Balloon	12-13.5-15	8	180	2.8	2/Box
FW24055	Rapide Esophageal Dilatation Balloon	15-16.5-18	8	180	2.8	2/Box
FW24065	Rapide Esophageal Dilatation Balloon	18-19-20	8	180	2.8	2/Box

INFLATION DEVICE STERILE

Order Number	Description	Package Units
DB25011	Micro-Tech Balloon Inflation Device 60 ml with Lock	1/Box
DB25012	Micro-Tech Balloon Inflation Device 60 ml with Lock	5/Case

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTEndoscopy.com
customercare@mtendoscopy.com

MPW23000 Rev.5