

ILLUMINATE | MAGNIFY | IDENTIFY PROTRAP LUMA



Low light environments can make identifying diminutive polyps a difficult task. The ProTrap Luma™ In-Line polyp trap comes with a built in light source and 2x magnification to help visualization.

The quick-change trays allow for fast and easy tray swaps to minimize suction downtime. The trays feature a bi-directional design which allows them to be flipped over to assist in debris clearing. When the sample is collected, the tray with sample can be directly placed in the specimen jar to minimize risks of tissue damage and biohazard exposure.

- Identify samples with confidence
- Reduce time off suction
- Prevent unnecessary sample handling



SPECIFICATIONS

PROTRAP LUMA [Ⓢ]

Order Number	Description	Package Units
TP54411	ProTrap Luma In-Line Polyp Trap	25/Box
TP54421	ProTrap Luma In-Line Polyp Trap Specimen Tray	50/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.