


# TARGETED THERMAL HEMOSTASIS ENSURE COAGULATION FORCEPS




EnSure™ reliable hemostasis of bleeding vessels without affecting surrounding tissue.

The EnSure Coagulation Forceps are designed to cauterize and coagulate GI bleeds. These single-use hemostatic forceps deliver targeted monopolar coagulation to isolate vessels and achieve hemostasis. Bleeding vessels can be approached with precision due to the reliable 360° bidirectional rotation.



Markings on catheter to help determine insertion length within the endoscope.



360° bidirectional rotation for targeting bleeds.

## SPECIFICATIONS

### ENSURE COAGULATION FORCEPS STERILE ②

Order Number	Description	Jaw Style	Working Length (cm)	Outer Diameter (mm)	Jaw Opening (mm)	Package Units
CF48361	Ensure Coag Forceps Grooved Jaw	Grooved Jaw	230	2.3	6.3	10/Box
CF48331	Ensure Coag Forceps Grooved Jaw	Grooved Jaw	165	2.3	6.3	10/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.