

## IDEAL BITE BIOPSY FORCEPS

With our wide selection of cost-effective disposable biopsy forceps, you are perfectly equipped for every procedure.

Available jaw diameters of 1.8 mm, 2.3 mm and 3.0 mm, and are offered in Oval, Serrated, Alligator jaws, with or without a spike

- Available lengths of 160 cm, 180 cm and 230 mm
- All forceps feature entry and exit indicator markings




HINGE



PRECISION JAWS

SINGLE DRIVEWIRE



CUSHIONED RINGS

THREE-RING HANDLE

SHORT THROW

### KEY BENEFITS

#### ERGONOMIC HANDLE

Three-ring comfort with cushioned lining and short throw provides improved grip and accommodates smaller hands.

#### DIRECT DRIVE SYSTEM

The single drive wire and hinge linkage deliver smooth activation around curves and tortuous anatomy. The system delivers more closure force to the jaws which provides for improved cutting ability and assures complete jaw closure.

#### SCOPE FRIENDLY DESIGN

The design of the Micro-Tech forceps family is based upon the robust and proven design of original reusable biopsy instruments. These designs provide optimal bite, while minimizing the risk of scope damage.

## BIOPSY FORCEP FOR 3.2 MM SCOPE CHANNELS STERILE ②

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40061	129-0177	Jumbo Biopsy	Oval	Blue	230	3.2	10/Box
BF40071	129-0178	Jumbo Biopsy with Spike	Oval	Blue	230	3.2	10/Box
BF40171	129-0182	Jumbo Biopsy	Alligator	Blue	230	3.2	10/Box
BF40181	129-0183	Jumbo Biopsy with Spike	Alligator	Blue	230	3.2	10/Box

## BIOPSY FORCEP FOR 2.8 MM SCOPE CHANNELS STERILE ②

Order Number	Henry Schein Item Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40081	129-0179	Standard Forceps	Oval	Orange	180	2.8	10/Box
BF40091	129-0677	Standard Forceps	Oval	Blue	230	2.8	10/Box
BF40096	128-5668	Standard Forceps	Oval	Blue	230	2.8	50/Box
BF40111	129-0180	Standard Forceps with Spike	Oval	Orange	180	2.8	10/Box
BF40121	129-0678	Standard Forceps with Spike	Oval	Blue	230	2.8	10/Box
BF40126	128-5669	Standard Forceps with Spike	Oval	Blue	230	2.8	50/Box
BF40191	129-0184	Standard Forceps	Alligator	Orange	180	2.8	10/Box
BF40201	129-0679	Standard Forceps	Alligator	Blue	230	2.8	10/Box
BF40206	128-5670	Standard Forceps	Alligator	Blue	230	2.8	50/Box
BF40221	129-0185	Standard Forceps with Spike	Alligator	Orange	180	2.8	10/Box
BF40231	129-0680	Standard Forceps with Spike	Alligator	Blue	230	2.8	10/Box
BF40236	128-5671	Standard Forceps with Spike	Alligator	Blue	230	2.8	50/Box
BF40431	132-9020	Standard Forceps	Serrated	Orange	180	2.8	10/Box
BF40441	132-9017	Standard Forceps	Serrated	Blue	230	2.8	10/Box
BF40446	131-2212	Standard Forceps	Serrated	Blue	230	2.8	50/Box
BF40411	132-9018	Standard Forceps with Spike	Serrated	Orange	180	2.8	10/Box
BF40421	132-9019	Standard Forceps with Spike	Serrated	Blue	230	2.8	10/Box
BF40426	131-2211	Standard Forceps with Spike	Serrated	Blue	230	2.8	50/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
customer@mtendoscopy.com