

## GRASP, APPROXIMATE, AND CLOSE DUAL ACTION TISSUE CLOSURE DEVICE



Gastrointestinal defect closure devices should be quick, easy to use, and effective.

The DAT™ Closure Device has a center post flanked by two arms which can operate independently to capture tissue and support manipulation to close mucosal defects.

The DAT closure device can be repositioned before deployment and can be used alone or to approximate tissue margins and hold them in place while other through-the-scope clips are placed to further secure defect closure.



### SPECIFICATIONS

#### DAT CLOSURE DEVICE STERILE

Order Number	Description	Opening Width (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
TC35511	DAT Closure Device - Gastroscope	15	165	3.2	Each
TC35521	DAT Closure Device - Colonoscope	15	195	3.2	Each
TC35531	DAT Closure Device - Colonoscope	15	235	3.2	Each

#### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.