

EFFICIENT RESECTION COLD SNARES

OVAL OR DIAMOND LOOP

- Designed specifically to resect diminutive polyps
- Available in two sizes – 10 mm and 15 mm
- Thin cutting wire provides for a guillotine-type cut
- Loop shape assists in capturing smaller polyps
- Stiffer catheter aids in insertion and provides support during resection

0.23 MM WIRE

HANDLE MAY
BE ROTATED

SMOOTH GLIDE

FIRM CATHETER

KEY BENEFITS

DEDICATED COLD SNARE

Specifically designed for resecting diminutive polyps using a cold snare technique, the oval or diamond shape loop and flexible wire capture smaller polyps with ease.

THIN WIRE

The thin 0.23 mm cutting wire provides a guillotine cut that reduces the force needed to complete the resection and may limit tearing or snagging.

COMFORTABLE ERGONOMIC HANDLES

Accommodating various hand sizes, the ergonomic three-ring handle features a shortened throw and comfortable grip. The stiff catheter provides support for a smooth, effortless cut.

SPECIFICATIONS

COLD SNARES STERILE 2

Order Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
CS50011	Diamond Thin Wire	10	7	230	2.8	10/Box
CS50021	Diamond Thin Wire	15	7	230	2.8	10/Box
CS50031	Oval Thin Wire	10	7	230	2.8	10/Box
CS50041	Oval Thin Wire	15	7	230	2.8	10/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

[MTEndoscopy.com](https://www.MTEndoscopy.com)
customer care@mtendoscopy.com

MPW50000 Rev.4