

ELECTROSURGICAL PRECISION GOLDKNIFE



The Micro-Tech Electrosurgical Knives have gold plating on the exposed metal posts, which may reduce tissue residue from adhering to the knife because gold is a highly effective conductor.

T KNIFE

The T tip GoldKnife has integrated electrosurgical and water injection functions, eliminating the need for instrument exchange which may save you time.

IT KNIFE

The IT tip GoldKnife features a ceramic tip which provides insulation, minimizing invasiveness while providing a 4.0 mm cutting electrode.

SPECIFICATIONS

GOLDKNIFE T STERILE

Order Number	Description	Working Channel (mm)	Length (cm)	Cutting knife length when extended (mm)
MK55131	GoldKnife T, Irrigating	2.8	195	1.5
MK55141	GoldKnife T, Irrigating	2.8	235	1.5
MK55151	GoldKnife T, Irrigating	2.8	195	2.0
MK55161	GoldKnife T, Irrigating	2.8	235	2.0
MK55171	GoldKnife T, Irrigating	2.8	195	4.0
MK55181	GoldKnife T, Irrigating	2.8	235	4.0
MK55211	GoldKnife IT	2.8	195	4.0
MK55221	GoldKnife IT	2.8	235	4.0

ACTIVE CORD

Order Number	Description	Package Units
QB51011	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.