

VISUALIZED DISECTION SURETRAC TRACTION SYSTEM



SureTrac™ retracts tissue to assist in tissue exposure, resection, and removal within the stomach and colon. The SureTrac system features a preloaded elastic silicone band on a primary clip and a secondary clip to anchor on tissue. The superelastic silicone band is designed to retract and release elastic potential energy for dynamic traction; this provides enhanced visualization and enables tissue management while dissecting submucosal lesions.

Each kit includes:

- A Primary Device, preloaded with a silicone band
- A Secondary Device with the same size clip as the Primary Device
- Available in either SureClip® or Lockado® clip design

SPECIFICATIONS

SURETRAC KIT STERILE

Order Number	Description	Clip Size	Clipping Device Style	Minimum Channel Size	Package Units
TS63041	SureTrac Elastic Traction Device	11 mm	SureClip	2.8 mm	1/Box
TS63061	SureTrac Elastic Traction Device	16 mm	SureClip	2.8 mm	1/Box
TS63101	SureTrac Elastic Traction Device	11 mm	Lockado	2.8 mm	1/Box
TS63121	SureTrac Elastic Traction Device	16 mm	Lockado	2.8 mm	1/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.