

PROTECT, COLLECT, CLEAN SINGLE-USE VALVE KITS



QL71121

Disposable air/water valves are intended to control inflow of gases and water and prevent backflow. Disposable suction valves are intended to enable suction and prevent the inflow of air. Single-use valves provide an alternative to manual cleaning and reprocessing of reusable valves, which can help to improve infection prevention practices.

- Compatible with Olympus & Pentax scopes
- Offers consistent performance by using a new valve kit for each procedure



QL71111, QL71112



QL71101

SPECIFICATIONS

SINGLE-USE VALVE KITS STERILE

Order Number	Description	Package Units
QL71151	Air/Water & Suction, Biopsy Valve Kit, Olympus, 3 Piece Kit	50/Box
QL71101	Air/Water & Suction, Biopsy Valve, Aux Water Kit, Olympus, 4 Piece Kit,	50/Box
QL71111	Air/Water & Suction, Biopsy Valve, Aux Water, Cleaning Adapter Kit, Olympus, 5 Piece Kit	10/Box
QL71112	Air/Water & Suction, Biopsy Valve, Aux Water, Cleaning Adapter Kit, Olympus, 5 Piece Kit, (4 Boxes of 10)	40/Case
QL71161	Air/Water & Suction, Biopsy Valve Kit, Pentax, 3 Piece Kit	50/Box
QL71121	Air/Water & Suction, Biopsy Valve, Aux Water Kit, Pentax, 4 Piece Kit	50/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.