


NEXT-GEN VISUALIZATION

EYEMAX SINGLE-USE CHOLANGIOSCOPE



Discover the advantages of enhanced resolution by using the eyeMAX™ Cholangioscope. With superior optical clarity and improved flexibility, eyeMAX provides unmatched access and clear visualization to aid diagnosis. Available in two sizes, the 3.9 mm O.D. scope includes a 2.0 mm working channel that accommodates our 1.6 mm eyeMAX biopsy forceps, potentially increasing biopsy yield. Meanwhile, the 3.2 mm version comes equipped with a 1.2 mm working channel that works seamlessly with the 1.0 mm eyeMAX biopsy forceps, enabling deeper access into the patient's anatomy.



CHOLANGIOSCOPES STERILE

Order Number	Version	Field of View	Insertion Diameter (mm)	Length (mm)	Working Channel Diameter (mm)	Package Units
CD63681	H	120°	3.9	220	2.0	1/Box
CD63701	C	120°	3.9	220	2.0	1/Box
CD63691	H	120°	3.2	220	1.2	1/Box
CD63761	C	120°	3.2	220	1.2	1/Box

DIGITAL CONTROLLER

Order Number	Voltage	Frequency	Functions
SW71212	120V	60Hz	LED, white color balance, brightness, zoom, image shape change

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.