

PRECISION RESECTION
**EDGEHOG
HOT SNARE**

EdgeHog™ is our most responsive hot snare with a flat monofilament wire which was designed to cut through tissue with precision and ease. It also features an ultra-stiff wire with a rounded silhouette and is rotatable, ideal for advanced EMR procedures.

REDUCED SLIPPAGE

The stiff monofilament wire indents into tissue and the flat wire sits perpendicular to the cutting plane which allows for maximum grip and hold on tissue. The proximal "V-shaped" design of the wire provides stability so the snare remains flat and can engage tissue. The result is less slippage as it closes on tissue.

FLAT MONOFILAMENT WIRE

The monofilament wire conducts energy efficiently. This helps the snare resect tissue efficiently, providing a clean margin with a razor-thin cautery effect which may reduce thermal damage of healthy tissue.



EDGEHOG HOT SNARE

Order Number	Description	Loop Diameter	Working Length	Package Units
PS51181	EdgeHog Hot Snare	20 mm	230 cm	10/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.