


NEXT-GEN VISUALIZATION  
**EYEMAX SINGLE-USE CHOLANGIOSCOPE**



Discover the advantages of enhanced resolution by using the eyeMAX™ Cholangioscope. With superior optical clarity and improved flexibility, eyeMAX provides unmatched access and clear visualization to aid diagnosis. Available in two sizes, the 3.9 mm O.D. scope includes a 2.0 mm working channel that accommodates our 1.6 mm eyeMAX biopsy forceps, potentially increasing biopsy yield. Meanwhile, the 3.2 mm version comes equipped with a 1.2 mm working channel that works seamlessly with the 1.0 mm eyeMAX biopsy forceps, enabling deeper access into the patient's anatomy.



**CHOLANGIOSCOPES** STERILE 

Order Number	Version	Field of View	Insertion Diameter (mm)	Length (cm)	Working Channel Diameter (mm)	Package Units
CD63681	H	120°	3.9	220	2.0	1/Box
CD63701	C	120°	3.9	220	2.0	1/Box
CD63691	H	120°	3.2	220	1.2	1/Box
CD63761	C	120°	3.2	220	1.2	1/Box

**DIGITAL CONTROLLER**

Order Number	Voltage	Frequency	Functions
SW71212	120V	60Hz	LED, white color balance, brightness, zoom, image shape change

**CAN'T FIND IT?**

Additional items may be available. Contact us if you can't find what you need.