


PRECISE CANNULATION

SPHINCTEROTOME



The Micro-Tech sphincterotome features a patent pending¹ opening mechanism to facilitate both short wire and long wire procedures with one device. The rotatable design and locking mechanism provide maximum control during cannulation and sphincterotomy.

KEY BENEFITS

FLEXIBILITY

Compatible with 0.035" guidewires, the rotatable tip can be actuated to ensure proper cannulation and sphincterotomy positioning.

CONTROL

The c-channel catheter allows for a splitting movement to separate wire from catheter for physician guided cannulation.

ATRAUMATIC TIP

The tapered atraumatic tip helps to ease insertion and tracking.

SHOW PRELOADED SPHINCTEROTOMES STERILE EO ⓧ


Order Number	Tip Length (mm)	Cutting Wire Insulation	Cutting Wire (mm)	Preloaded Wire Length (cm)	Working Length (cm)	Package Units
SS75041	5	Insulated	20	260	200	1/Box
SS75061	5	Insulated	20	450	200	1/Box
SS75101	5	Insulated	25	260	200	1/Box
SS75121	5	Insulated	25	450	200	1/Box*
SS75161	5	Insulated	30	260	200	1/Box
SS75181	5	Insulated	30	450	200	1/Box*

Preloaded with 0.035" Hydro-Slide Straight Tip Guidewire.

SHOW SPHINCTEROTOMES STERILE EO ⓧ

Order Number	Tip Length (mm)	Cutting Wire Insulation	Cutting Wire (mm)	Working Length (cm)	Package Units
SS75011	5	None	20	200	1/Box*
SS75021	5	Insulated	20	200	1/Box
SS75071	5	None	25	200	1/Box*
SS75081	5	Insulated	25	200	1/Box
SS75131	5	None	30	200	1/Box*
SS75141	5	Insulated	30	200	1/Box

*Non-Stock Item ¹US patent pending



MAXIMUM FLEXIBILITY

BILIARY STONE RETRIEVAL BALLOON

Micro-Tech's stone extraction balloon is comprised of a high-grade latex for reliability during endoscopic bile stone extraction. The 3 stage balloons come in above and below injection options to meet a wide range of clinical needs. The catheter features a 3-lumen design for insufflation, guidewire access, and contrast injection.

KEY BENEFITS

PROCEDURAL FLEXIBILITY

3-stage balloons are sized in 3 mm increments between 9 and 21 mm.

INSERTION AND TRACKING

The tapered atraumatic tip helps to ease insertion and improve tracking along the wire into biliary anatomy.

COMPATIBILITY

The balloon catheter is compatible with 0.035" guidewires.

STONE RETRIEVAL BALLOONS



Order Number	Catheter Diameter (Fr)	Balloon Stages	Injection Port Location	Working Length (cm)	Package Units
SB75921	7	9/12/15 mm	Above (distal to balloon)	200	1/Box
SB75931	7	12/15/18 mm	Above (distal to balloon)	200	1/Box
SB75941	7	15/18/21 mm	Above (distal to balloon)	200	1/Box
SB76131	7	9/12/15 mm	Below (proximal to balloon)	200	1/Box*
SB76141	7	12/15/18 mm	Below (proximal to balloon)	200	1/Box*
SB76151	7	15/18/21 mm	Below (proximal to balloon)	200	1/Box*

*Non-Stock Item

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.

