


ACCELERATED ACCESS  
**HYDRO-SLIDE  
GUIDEWIRE**

A reliable guidewire is the core of a successful procedure. The Hydro-Slide Guidewire provides the stiffness, visibility, and atraumatic tip needed to bridge strictures in the GI tract.

The PTFE jacket, Hydrophilic, and Nitinol core provide improved pushability while the soft tungsten loaded tip allows for high visibility under fluoroscopy. The yellow and blue striped jacket is an easy-to-see indicator when advancing the guidewire under direct visualization.

- Straight and Angled tips
- 260 cm and 450 cm working lengths
- Individually packaged, 2 per box



## HYDRO-SLIDE .025" GUIDEWIRE STERILE X-RAY

Order Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10011	Hydro-Slide Guidewire Straight	.025	260	2/Box
GW10021	Hydro-Slide Guidewire Angled	.025	260	2/Box
GW10031	Hydro-Slide Guidewire Straight	.025	450	2/Box*
GW10041	Hydro-Slide Guidewire Angled	.025	450	2/Box*

## HYDRO-SLIDE .035" GUIDEWIRE STERILE X-RAY

Order Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10051	Hydro-Slide Guidewire Straight	.035	260	2/Box
GW10071	Hydro-Slide Guidewire Angled	.035	260	2/Box
GW10091	Hydro-Slide Guidewire Straight	.035	450	2/Box*
GW10111	Hydro-Slide Guidewire Angled	.035	450	2/Box*

## HYDRO-SLIDE .035" STIFF GUIDEWIRE STERILE X-RAY

Order Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10061	Hydro-Slide Guidewire Stiff Straight	.035	260	2/Box
GW10081	Hydro-Slide Guidewire Stiff Angled	.035	260	2/Box
GW10101	Hydro-Slide Guidewire Stiff Straight	.035	450	2/Box*
GW10121	Hydro-Slide Guidewire Stiff Angled	.035	450	2/Box*

\*Non-Stock Item

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.