

SHOW:

# GUIDE WIRE LOCKING DEVICE

The SHOW Short Wire Locking Device simplifies and enhances the process of guidewire management during ERCP procedures. By fixing the guidewire in place, the SHOW Short Wire Locking Device allows for seamless device exchanges while minimizing the risk of wire displacement.



Featuring a **One-click** mounting system allowing for quick and easy attachment to the working channel



Can hold up to **two guidewires in parallel**, offering greater flexibility during device exchanges



**Double membrane** valve aids in the reduction of air and bile reflux, to promote a cleaner and more controlled procedure environment.

## GUIDE WIRE LOCKING DEVICE

Order Number	Description	Package Units
GL62781	SHOW Guidewire Locking Device, Oly, Fuj	10/Box
GL62791	SHOW Guidewire Locking Device, Pentax	10/Box

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.