

ELECTROSURGICAL PRECISION

GOLDKNIFE



The Micro-Tech Electrosurgical Knives have gold plating on the exposed metal posts, which may reduce tissue residue from adhering to the knife because gold is a highly effective conductor.

The T-tip GoldKnife™ has integrated electrosurgical and water injection functions, eliminating the need for instrument exchange which may save you time.



GOLDKNIFE T STERILE

Order Number	Description	Working Channel (mm)	Length (cm)	Cutting knife length when extended (mm)	Package Units
MK55131	GoldKnife T, Irrigating	2.8	195	1.5	1/Box
MK55141	GoldKnife T, Irrigating	2.8	235	1.5	1/Box
MK55151	GoldKnife T, Irrigating	2.8	195	2.0	1/Box
MK55161	GoldKnife T, Irrigating	2.8	235	2.0	1/Box

ACTIVE CORD

Order Number	Description	Package Units
QB51011	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box