

ENHANCED SAFETY AND EFFICIENCY **SINGLE-USE ENDOSCOPE VALVES AND ACCESSORIES**



Reusable medical devices, like endoscope valves, can present numerous challenges for patients and clinicians, as they can become contaminated, posing an infection risk if not adequately reprocessed between patient uses. Additionally, the manual cleaning and reprocessing of these devices can be time-consuming, involving multiple steps to ensure complete sterilization and traceability to a specific endoscope.

Designed to address these challenges, our disposable valves and accessories are a reliable alternative to reusable devices that can help streamline procedures and enhance infection prevention practices by eliminating the need for timely cleaning, reprocessing, and traceability procedures and reduce the potential risk of infection from cross-contamination.

COMPREHENSIVE OFFERINGS TO MEET YOUR NEEDS

Our single-use valves and accessories are compatible with Olympus®, Fujifilm®, and Pentax® Endoscopes, ensuring seamless integration with your existing equipment. Available in 3-, 4-, or 5-piece kits and single-pack configurations, our comprehensive offerings feature air/water valves, suction and biopsy valves, auxiliary water connectors, and cleaning adapters, providing the right solution to meet your needs.

OLYMPUS® ENDOSCOPES



FUJIFILM® 700 SERIES ENDOSCOPES



FUJIFILM® 500/600 SERIES ENDOSCOPES



PENTAX® ENDOSCOPES



SINGLE-USE ENDOSCOPE VALVES & ACCESSORIES STERILE

| Order Number | Description | Package Units |
|--|---|---------------|
| OLYMPUS® ENDOSCOPES | | |
| QL72111 | 5-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector & Cleaning Adapter) | 50/Box |
| QL72101 | 4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector) | 50/Box |
| QL72151 | 3-Piece Kit (Air/Water, Suction, & Biopsy Valves) | 50/Box |
| QL72191 | Single-Pack Biopsy Valve (Olympus®/Fujifilm® Compatible) | 100/Box |
| QL72171 | Single-Pack Auxiliary Water Connector | 100/Box |
| QL72231 | Single-Pack Cleaning Adapter* | 100/Box |
| FUJIFILM® 700 ENDOSCOPES | | |
| QL72131 | 5-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector & Cleaning Adapter) | 50/Box |
| QL72281 | 4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector) | 50/Box |
| QL72301 | 3-Piece Kit (Air/Water, Suction, & Biopsy Valves) | 50/Box |
| QL72191 | Single-Pack Biopsy Valve (Olympus®/Fujifilm® Compatible) | 100/Box |
| QL72211 | Single-Pack Auxiliary Water Connector | 100/Box |
| QL72261 | Single-Pack Cleaning Adapter* | 100/Box |
| FUJIFILM® 500/600 SERIES ENDOSCOPES | | |
| QL72291 | 5-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector & Cleaning Adapter) | 50/Box |
| QL72141 | 4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector) | 50/Box |
| QL72271 | 3-Piece Kit (Air/Water, Suction, & Biopsy Valves) | 50/Box |
| QL72191 | Single-Pack Biopsy Valve (Olympus®/Fujifilm® Compatible) | 100/Box |
| QL72211 | Single-Pack Auxiliary Water Connector | 100/Box |
| QL72251 | Single-Pack Cleaning Adapter* | 100/Box |
| PENTAX® ENDOSCOPES | | |
| QL72121 | 4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector) | 50/Box |
| QL72161 | 3-Piece Kit (Air/Water, Suction, & Biopsy Valves) | 50/Box |
| QL72201 | Single-Pack Biopsy Valve | 100/Box |
| QL72181 | Single-Pack Auxiliary Water Connector | 100/Box |

*Single-pack cleaning adapter is non-sterile.

All trademarks are the property of their respective owners.

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.